


# Huygens Essential User Guide



*Scientific Volume Imaging*





# ***Huygens Essential User Guide***

Scientific Volume Imaging b.v.

Alexanderlaan, 14  
1213 XS Hilversum

P.O. box 615  
1200 AP Hilversum

The Netherlands

<http://www.svi.nl>

Copyright © 1995–2006 by Scientific Volume Imaging b.v.  
Alexanderlaan, 14. 1213 XS Hilversum  
P.O. Box 615. 1200 AP Hilversum  
The Netherlands

All rights reserved

Cover illustration:


Macrophage recorded by Dr. James Evans (White-head Institute, MIT, Boston MA, USA) using widefield microscopy, as deconvolved with Huygens.

At the right the same dataset again: macrophage fluorescently stained for tubulin (yellow/green), actin (red) and the nucleus (DAPI, blue). Left part: original data; right part: as deconvolved with the classical Maximum Likelihood Estimation method (MLE).

The image was visualized using the Simulated Fluorescence Process (SFP) volume rendering package from Scientific Volume Imaging.

Image in figures 22, 24, 27 and 29: isolated Rat Hepatocyte couplet recorded by Dr. Permsin Marbet at the Department of Anatomy of the University of Basel, Switzerland (head: Prof. Lukas Landmann), as deconvolved with Huygens.

Image in figure 28: FISH-stained cell nucleus recorded at the Nuclear Organization Group, SILS, University of Amsterdam, The Netherlands (head: Prof. Roel van Driel).



This document was typeset with [OpenOffice.org](http://OpenOffice.org).

# Table of Contents

<b>1. What is Huygens Essential?</b>	<b>1</b>
<b>2. Installing Huygens Essential</b>	<b>1</b>
Mac OS X	1
Microsoft Windows	1
Linux	2
SGI Irix 6.5	2
IBM AIX 5.2	2
After the installation	2
The license string	2
Obtaining a license string	3
Installing the license string	3
Location of the license file	3
Trouble shooting license strings	4
Updating the software	4
Removing the software	4
System requirements for Huygens Essential	5
Windows and Linux	5
Mac OS X	5
SGI Irix, IBM AIX	5
Running the program	5
Support on installation	6
<b>3. The image restoration process</b>	<b>6</b>
The processing stages	6
Loading an image	7
Preprocessing	8
Converting a data set (optional)	8
Time series	8
Adapting the image	8
Processing brightfield images	8
Setting the image channel colors	9
Image statistics	9
Verifying the microscopic parameters	9
Parameters Templates	10
The intelligent cropper	11
Cropping an image in x y z	11
Cropping an image in time	12
Removing channels	12
Stage 1: Parameter tuning	12
Stage 2: The image histogram	12
Stage 3: Estimate the average background in the image	13
Stage 4: Deconvolution	14
The Point Spread Function (PSF)	14
The Classical Maximum Likelihood Estimation (CMLE) algorithm	15
Number of iterations	15
Signal to noise ratio	15
Quality threshold	15
Iteration mode	15
Bleaching correction	15
Stopping the MLE algorithm	16
Finishing or restarting a deconvolution run	16
z-drift corrector for time series	16
Saving the result	16
Saving the restored image	16
Saving your Task report	17
Multi channel images	17
Deconvolving a channel in a multi channel image	17
Joining the results	17
Batch processor	18

Start up	18
Window description	18
Usage	19
Adding one task	19
Adding multiple tasks	19
Running the batch job	20
Menus	20
File menu	20
Task menu	20
More information	20
<b>4. Huygens Essential Visualization Tools</b>	<b>21</b>
<i>The Twin Slicer</i>	21
Color mode	21
Contrast mapping mode	22
Time series	22
<i>The MIP Renderer</i>	23
Soft threshold	23
Rendering a movie	23
<i>The SFP Renderer</i>	24
Summary	24
SFP fundamentals	24
Rendering a movie	26
<i>The Surface Renderer</i>	26
Hue Selector	27
Transparency depth	27
<b>5. Huygens Essential Analysis Tools</b>	<b>28</b>
<i>The Object Analyzer</i>	28
How to use the Object Analyzer	28
Mouse mode	29
Frame selector	29
Render pipes	30
Object removal	30
Regions of interest (ROI)	30
Data analysis	31
Visualization parameters	32
Save the results	32
Read more	32
<i>The Colocalization Analyzer</i>	32
Iso-colocalization object analysis	33
How to use the Colocalization Analyzer	33
Backgrounds vs. thresholds in colocalization	35
Read more	35
<b>6. The PSF Distiller</b>	<b>36</b>
Beads for PSF distillation	36
<i>The Distiller stages</i>	37
Starting the Distiller	37
Parameter check stage	37
Averaging stage	37
Confocal and two photon bead images	37
Widefield bead images	38
Distillation stage	38
Assembly stage	38
<b>7. Establishing image parameters</b>	<b>38</b>
<i>Image size</i>	38
Brick wise processing	39
<i>Signal to Noise Ratio (SNR)</i>	39
<i>Blacklevel</i>	40
<i>Sampling densities</i>	40
<i>Computing the backprojected pinhole radius</i>	41
Airy disk as unit for the backprojected pinhole	42
Converting from integer parameter	42
Shape factor	42
Airy disk as unit for the backprojected radius of a square pinhole	43

Computing the backprojected pinhole distance in Nipkow spinning disks	43
<i>Pinhole radius tables</i>	43
Leica confocal microscopes	43
TCS 4d, SP1, NT	43
TCS-SP2	44
Zeiss confocal microscopes	44
Olympus confocal microscopes	45
Biorad confocal microscopes	45
Checking the Biorad system magnification	45
A supplied calibration curve	46
<b>8. Improving the quality of your images</b>	<b>46</b>
<i>Data acquisition pitfalls</i>	46
Refractive index mismatch	46
Clipping	46
Undersampling	47
Do not undersample to limit photodamage	47
Bleaching	47
Illumination instability	47
Mechanical instability	47
Thermal effects	48
Internal reflection	48
<i>Deconvolution improvements</i>	48
Acquire an experimental PSF	48
Spherical aberration correction	48
Improve the estimated parameters	48
<b>9. Appendix</b>	<b>49</b>
<i>License string details</i>	49
<i>Questions</i>	50
Where can I find support on the web?	50
What does the quality factor mean while running Huygens?	50
Can I deconvolve a TIFF series?	51
TIFF file series naming convention	51
Can I deconvolve a single plane widefield image?	52
Can I deconvolve a single TIFF image?	52
How do I generate a debug log?	52
<i>Addresses and URLs</i>	53
Where can we be reached?	53
Scientific Volume Imaging b.v.	53
Distributors	53
Support and FAQ	53
Knowledge database (FAQ)	53
The SVI-wiki	53
Starting points	53
Quick reference	54
<b>Alphabetical Index</b>	<b>55</b>





## 1. What is Huygens Essential?

Huygens Essential is an image processing software package tailored for restoration, visualization and analysis of microscopic images. Its wizard driven user interface guides you through the process of deconvolving images from light microscopes. Huygens Essential is able to deconvolve a wide variety of images ranging from 2D widefield (WF) images to 4D multi-channel multi-photon confocal images. To facilitate comparison of raw and deconvolved data or results from different deconvolution runs the Essential is equipped with a dual 4D slicer tool. You can also render 3D images and animations with its powerful visualization tools. Post-restoration analysis is possible using the interactive analysis tools.

Based on the same image processing engine (the compute engine) as Huygens Professional, Huygens Essential combines the quality and speed of the algorithms available in Huygens Professional with the ease of use of a wizard driven intelligent user interface.

Huygens Essential is based on cross-platform technology. It is available on various Microsoft Windows operating systems, Linux for Intel or AMD based systems, MAC OS X, IBM AIX and SGI's Irix 6.5. For AIX and Irix 64 bit multiprocessing versions are also available.

## 2. Installing Huygens Essential

You can download Huygens Essential from the SVI website: <http://www.svi.nl>


### Mac OS X

Go to the folder where you downloaded the distribution and double-click on it. It will be extracted by "StuffIt Expander" to a .pkg file, which will be placed in the same directory. Double click this .pkg file, and follow the installation wizard.

### Microsoft Windows

You have received an executable file for installation, for instance a file named Huygens280.exe. Place this file on your desktop and double click its icon to start the installation. During installation the directory

C:\Program files\SVI\ will be created by default. After completion the Huygens Essential and Huygens Scripting icons appear on your desktop. Double clicking on the Huygens Essential icon starts the program; the start-up window will be displayed (Figure 1).



**Figure 1. The start-up window on Microsoft Windows.** If no license string is installed the software runs in 'Freeware mode'. You can find your computer's ID number by using **HELP > ABOUT**

## **Linux**

The Huygens Essential Linux distribution is a 'rpm' file, for instance `huygens-2.7.0-p7.rpm`. Open a Unix shell, go to the directory where this file is located, become superuser and type: `rpm -ivh --force huygens-2.7.0-p7.rpm`

After installing the software type `essential` in a shell to start the software. A directory `/usr/local/svi` will be created; initialization scripts will be installed in `/usr/local/bin`.

## **SGI Irix 6.5**

Currently the Irix distribution is a single 'tardist' file containing various components. By default all components are installed. Become superuser and type:

```
swmgr -f dist65-2.7.0-p7.tardist
```

Press **START** in the Software Manager window.

After installing the software type `essential` in a shell to start the software. The program will display the start-up window (Figure 1). A directory `/usr/local/svi` will be created; the executables will be installed in `/usr/local/bin`.

## **IBM AIX 5.2**

Log in as root on your workstation or ask your system administrator to do so. Huygens for AIX is distributed in a tar file. Go to the directory where you downloaded it to. After unpacking the distribution file with `tar xvf my_aixfile.tar` you will find three new files: `svi.tar.Z`, `README.AIX` and a shell script `AIX-install.sh`. To install the software you have to run the file `AIX-install.sh` by typing `./AIX-install.sh`. This script will unpack the `svi.tar.Z` file under `/usr/local` and then will do some post installation tasks and verification. After a successful installation it will print the message "OK." to the screen.

A directory `/usr/local/svi` will be created; the executables will be installed in `/usr/local/bin`.

## **After the installation**

After a first-time installation there is not yet a license available. Still, you can start the software. Without a license it will run in 'Freeware mode'. Among others, this gives you access to the License tools in the **HELP** menu. The next section explains how to obtain and install a license string.

On AIX, Irix and Linux start the software by clicking its icon or by typing `essential` into a Unix shell. On Mac OS X and Windows click its icon. The software will open in Freeware mode and display the start-up window (Figure 1).

## **The license string**

The license key used by all SVI software is a single string per licensed package. It may look as follows:

```
HuEss-2.7-wcnp-d-tv-emnps-eom2008Dec31-e7b7c623393d708e-  
{frank@svi.nl}-4fce0dbe86e8ca4344dd
```

At startup Huygens Essential searches for a license file `huygensLicense` which contains a license string. This license string is provided by SVI via e-mail. Installing the license string is the same for all platforms, though on Linux, Irix and AIX only the superuser can do this.

## Obtaining a license string

If you are not upgrading from a previous installation it is likely that a license is not yet available. To enable us to generate a license string for you we need the *fingerprint* of your computer, the *system ID* number. If you have not done so already, start Huygens Essential. The system ID is displayed in the **HELP > ABOUT** dialog (Figure 2). Send it to [sales@svi.nl](mailto:sales@svi.nl), and you will be provided with a license string. To prevent any typing error use the **COPY** button to save the ID to the clipboard. You can print it into your mail message with the **EDIT > PASTE** menu item of your mail program.

In this dialog box you can also find a button to **CHECK FOR HUYGENS UPDATES** on our company server.

## Installing the license string

Select the license string in your email message and copy it to the clipboard using **EDIT > COPY** in your mailing program. Start Huygens Essential and go to **HELP > LICENSE**: a dialog box pops up (Figure 3). Then press the **ADD A NEW LICENSE** button and a new window will pop up (Figure 4). Paste the string into the text field using your keyboard<sup>1</sup>. Complete the procedure using the **ADD LICENSE** to add the string to the `huygensLicense` file. Please try to avoid typing the license string by hand: any small typing error will invalidate the license. With an invalid license, the software will remain in Freeware mode.




Figure 2 The **HELP > ABOUT** window. The system ID is shown at the bottom.




Figure 3. The **License manager dialog**. Left: as displayed on Mac OS X; right: as displayed on Irix. The License manager allows you to add, delete and troubleshoot licenses.




Figure 4. The **Add License dialog box**.

## Location of the license file

The license string is added to the file `huygensLicense` in the `svi` directory. On the different supported platforms this is located in:

- 1 Use your operating system's generic copy / paste operations: *Mac OS X*: `apple-c` / `apple-v`; *Windows*: `control-c` / `control-v`; *Linux, Irix, Aix*: on most common desktops copying is done simply by marking the text area with the mouse, and pasting by either middle mouse click or `control-v`.

- AIX, Irix and Linux: /usr/local/svi
- Mac OS X: depends on where you installed the software. A typical example is /Applications/SVI
- Windows: C:\Program Files\SVI\

On AIX, Irix and Linux and Mac OS X an alternative location is the user's home directory. On OS X this is especially convenient when updating frequently.



Restart Huygens Essential to activate the new license. This will enable the deconvolution or PSF distiller functionality.

### Trouble shooting license strings

The license string as used by SVI has the same appearance on all supported platforms. For each product<sup>2</sup> you need to have a license string installed. Select a license string in the License manager and press the EXPLAIN LICENSE button. All details for the current license will be listed. If you run into licensing problems you may use this information to analyze the problem. See *License string details* on page 49.




Figure 5. *The License details.*

### Updating the software

Download the new version from the SVI website at <http://www.svi.nl>. Proceed with the installation as explained above.



Do not uninstall the old version as this will delete your license string. On Mac OS X make sure you make a backup of the license string in a safe place before you remove the previous installation.

### Removing the software

- Irix: Open the Software Manager, select MANAGE INSTALLED SOFTWARE, and mark the packages you wish to remove.
- Linux: To remove Huygens Essential as root: `rpm -e huygens`
- Mac OS X: Drag the installation to the waste basket.
- Microsoft Windows: Clicking START in your Windows desktop and select: PROGRAMS > HUYGENS ESSENTIAL > UNINSTALL.



Removing the software will also cause your license string to be removed. If you prefer to uninstall your current version prior to installing a newer one, be sure to store your license string in a safe place.

<sup>2</sup> Huygens Essential, Huygens Professional and Huygens Scripting.

## System requirements for Huygens Essential

### Windows and Linux

Huygens Essential and Huygens Scripting run on the following Windows operation systems: Windows 2000, NT, 2003 Server and XP.

Linux: RedHat and SuSE distributions. Since Linux versions evolve rapidly best consult SVI's <http://www.svi.nl> web page to see which Linux distributions are currently supported.

A standard Ethernet card is required to provide your computer with a system ID.

- Processor: Pentium III or IV (Intel) or Athlon (AMD).
- Recommended RAM memory: 512 MB or larger (to run larger images like 512\*768\*50 voxels).
- Graphics card: any fairly modern card will do.

### Mac OS X

Mac OS X 10.2 or higher running on a G4 processor or higher. 512 MB or more RAM.

### SGI Irix, IBM AIX

Huygens Essential and Huygens Scripting run on all SGI equipment running Irix 6.5 on a MIPS R5000 processor or higher, IBM AIX 5.2 equipment with a Power4 processor or higher. The recommended RAM size is 512 MB or larger.

---

## Running the program

### Mac OS X

You can find the program icon in the installation directory. Clicking it will start Huygens Essential and open the main window (Figure 6). You can also run the program by typing `essential` at a shell prompt.

### Windows

The Windows installation procedure has automatically placed an Huygens Essential icon on your desktop. Clicking it will start Huygens Essential and open the main window (Figure 6).

### Linux

On Linux you can start Huygens Essential by typing `essential` at a shell prompt. It will start Huygens Essential and opens the main window (Figure 6).

If the shell is unable to find this command then typing the full path should help:

```
/usr/local/bin/essential
```

If this still does not help then Huygens Essential has not been installed correctly.

In Linux KDE desktop you may also start Huygens Essential from the Application menu and in the Gnome desktop from the Main Menu.

### Irix

On Irix you can start Huygens Essential by typing `essential` at a shell prompt which will start Huygens Essential and opens the main window (Figure 6).

If the shell is unable to find this command then typing the full path should help:

```
/usr/sbin/essential
```

If this still does not help then Huygens Essential has not been installed correctly.

### Adding to the shell search path



Users of the `csh` or `tcsh` shell can add the `/usr/local/bin` or `/usr/sbin/` directory to their shell search path by adding the following line to the `.cshrc` file in their home directory:

```
set path=(/usr/local/bin $path) or set path=(/usr/sbin $path)
```

You can inquire your shell by typing: `echo $shell`

### Support on installation

If you find any problem in installing the program or the licenses that you could not solve with the guidelines here included, please search the knowledge database or contact SVI on the addresses on page 53.




Figure 6. **The start-up window.** The gray area in the upper part is the work area where the thumbnail representation of the original image and its individual channels will be placed. The blue window in the left will show help text, task reports and supporting information during the various processing stages. In the gray bottom right field different dialog boxes will appear during processing. Also intermediate deconvolution results are shown here. The bottom bar is a status area.

## 3. The image restoration process

### The processing stages

Huygens Essential guides you through the process of microscopic image deconvolution (also referred to as 'restoration') in several stages. Each stage is composed of one or more tasks. While proceeding each stage is briefly described in the bottom-left *Task Info* window pane. The stages progress is indicated at the right side of the status bar (see Figure 15 on page 13). Additional information can be found in the `HELP > QUESTIONS` and `HELP > DICTIONARY` (Figure 7).

The following steps and stages are to be followed:

- Opening an image. A demo image (the `fab64.ics/ids` file pair) is placed in the `images` subdirectory of `svi`.
- stage P: The Preprocessing stage: loading an image, converting data sets, parameter check and cropping.
- stage 1: Tune the parameters. This stage will be skipped when you entered from the preprocessing stage. In the preprocessing stage you have already checked the parameter settings for the intelligent cropper - intelligent since it uses a-priori knowledge for setting the optimal cropping boundaries automatically. You will enter this stage from the latest one when clicking the <<RESTART button. This is useful if you wish to fine tune your parameters for the best deconvolution result, in particular when you like to set your parameters slightly different when using multi-channel images (see page 17).
- stage 2: Inspecting the image histogram.
- stage 3: Background estimation.
- stage 4: The deconvolution run.
- saving the result.




Figure 7. **The Dictionary.** The **DICTIONARY** and the **QUESTIONS** from **HELP** give additional information.

The different stages will be explained below

### Loading an image

Select **OPEN** from the **FILE** menu to enter the file browser and move to the directory where your images are stored. Select the image to be deconvolved, e.g. the `fab64.ics/ids` file pair in the `images` subdirectory of `svi`.

Several formats from microscope vendors are supported. If you have **TIFF** images to be processed please read *TIFF file series naming convention* on page 51 for the naming convention in order to be able to read a multi-dimensional image as a whole.

When the file is read successfully you can either press **START DECONVOLUTION** to begin processing your image or you can convert your data set with the **TOOLS** button.

If you have loaded a bead image you also can proceed selecting **START PSF DISTILLER** and proceed with generating a Point Spread Function (PSF, see page 14) from measured beads (See *The PSF Distiller* on page 36).



A special license is needed in order to launch the PSF Distiller.

You can **OPEN ADDITIONAL** images for reference purposes, but only the one named 'original' will be deconvolved during the guided restoration.

## Preprocessing

### Converting a data set (optional)

Before you press the `START DECONVOLUTION` button you can convert a 3D stack into time series images (Convert XYZ to XYZT) or vice versa, or you can convert a 3D stack into a time series of 2D images (XYZ to XYT) or vice versa. These functions can be found in the `TOOLS`.

Hint: If you have a data stack that is poorly sampled in *z* (not fulfilling the Nyquist criterion) you better interpret the different planes as independent (i.e. as 2D images) and do 2D deconvolution planewise while taking the optimal Nyquist criterion for *z* as imposed by the optical parameters (see the diagram in , page ). To do this in one run for all the planes, convert the 3D stack to a 2D-time series, do the deconvolution run, and convert back from 2D-time to 3D.

### Time series

A time series is a sequence of images recorded along time at uniform time intervals. Every recorded image is a time frame. The Huygens Essential is capable of automatic deconvolution of 2D-time or 3D-time data. There are some tools that are intended only for time series, as the confocal bleaching corrector or the *z*-drift corrector (page 16).

### Adapting the image

In the `TOOLS` menu you can find a contrast inverter helpful for the processing of brightfield images, see below.

A `CROP` tool is also available, but its use is recommended only after properly tuning the image parameters and will be explained in a later stage.



In the `TOOLS` menu you can also find a `MIRROR ALONG Z` tool, to flip the image when the coverslip is in the top. This is specially important in case of a refractive index mismatch: see *Spherical aberration correction* on page 48.

### Processing brightfield images

Brightfield imaging is not a 'linear imaging' process. In a linear imaging process the image formation can be described as the linear convolution of the object distribution and the point spread function (PSF, see page 14), hence the name deconvolution for the reverse process. So in principle one cannot apply deconvolution based on linear imaging to non linear imaging modes like brightfield and reflection. One could say the image formation in these cases IS linear because it is governed by linear superposition of amplitudes. However, microscopes do not measure light amplitudes but rather intensities, the absolute squared values of the amplitudes. Taking the absolute square destroys all phase information one would need to effectively apply deconvolution. Fortunately, in the brightfield case the detected light is to a significant degree incoherent. Because in that case there are few phase relations the image formation is largely governed by the addition of intensities, especially if one is dealing with a high contrast image.

In practice one goes about deconvolving brightfield images by inverting them (using `TOOLS > INVERT IMAGE`) and processing them further as incoherent fluorescence widefield images. Still, one should watch out sharply for interference like patterns (periodic rings and fringes around objects) in the measured image. As a rule these become pronounced in low contrast images. After the deconvolution run you may reverse to the original contrast setting.



### Setting the image channel colors

The Color picker tool (TOOLS > SET CHANNEL COLORS) allows you to alter the colors for the different image channels (see *Multi channel images* on page 17). The color of a particular channel can be edited by clicking the corresponding button. This opens a platform specific color editor (Figure 8).




Figure 8. **The Image Color Picker.** Clicking on a channel button (top) opens a color editor (right) by which the color can be changed. The color editor is platform specific. The picture on the right shows the Mac OS X color editor.

### Image statistics

Right-mouse-click on a thumbnail image and select SHOW PARAMETERS and you will find, besides the parameter settings, statistical information of the particular image. Amongst them are the mean, sum, standard deviation, norm, and position of the center of mass.

### Verifying the microscopic parameters

Next to the basic voxel data the Huygens Essential also tries to read as much as possible information about the microscopic recording conditions. However, depending on the file type this information may be incomplete or incorrect. In this first stage all parameters relevant for deconvolution are displayed and can be modified:

Optical parameters (first page):

- Microscope type
- Lens and medium refractive index
- Numerical Aperture (NA)

Optical parameters (second page):

- Backprojected pinhole radius in nm. 'Backprojected' means the size of the pinhole as it appears in the specimen plane, see *Computing the backprojected pinhole radius* on page 41.
- Backprojected distance between the pinholes in microns (only visible if the microscope type is 'Nipkow').
- Excitation and emission wavelengths
- Photon count (number of excitation photons involved in the fluorescence)
- Voxel sizes in the three directions  $x$ ,  $y$ ,  $z$  (third page)
- Summary of all parameters now in effect

If values are displayed in a **red background**, they are highly suspicious. An **orange background** indicates a non-optimal situation. Oversampling is also indicated with a **cyan background**, that becomes **violet** when it is very severe. Figure 9 shows settings which do not fulfill the criterion for the critical sampling distance versus numerical aperture. (See *Sampling densities* on page 40).

You can see and correct the image parameters not only at this deconvolution stage, but also at any time by right-clicking on the image thumbnail and selecting **SHOW PARAMETERS** or **CORRECT PARAMETERS**. Correcting will pop up the window shown on Figure 11.

### Parameters Templates

Once the proper parameters have been set and verified, they can be saved to a Huygens template file (suffix .hgst). Those template files are loadable in the very start of the wizard, hence the user can skip the parameters verification stage, provided that an image is to be restored with the same optical properties as the ones which were recorded on the template.




Figure 9. **Parameters check stage: Sampling.** Red coloring indicates a suspicious value, and orange a non-optimal value.




Figure 10. **Importing a microscopical parameters template.**

The **IMPORT MICROSCOPIC TEMPLATE** button will allow you to choose a template from a list of pre-saved template files which reside both in the common templates directory and in the user's personal template directory. The Huygens common templates directory is named **Templates**, and resides in the Huygens installation directory, namely `/usr/local/svi/Templates` on Unix systems, `C:\Program Files\SVI\Templates` on Windows and `/Applications/SVI/Templates` on the Mac OS X. The user's personal templates directory is called `.svi` on the Unix platforms and `SVI` on Windows, and it can be found in the user's home directory on Unix, and in `C:\Documents and Settings\user_name` on Windows. You can also choose to load a template file from a different location by pressing the **FROM DISK...** option.




Figure 11. **Microscopical parameters corrector.**

The EXPORT TEMPLATE button will allow you to either save the template to the personal template directory by choosing the TO DISK... option, or overwrite one of the existing templates by selecting them from the list.

The Huygens template is a simple .xml file which can be edited 'by hand' as well (see examples in the common Templates directory).

### The intelligent cropper

The time needed to deconvolve an image increases more than proportional with its volume. Therefore, deconvolution can be accelerated considerably by *cropping* the image. Huygens Essential is equipped with an intelligent cropper which automatically surveys the image to find a reasonable proposal for the crop region. In computing this initial proposal the microscopical parameters are taken into account, making sure that cropping will not have a negative impact on the deconvolution result. Because the survey depends on accurate microscopical parameters it is recommended to use the intelligent cropper as final step in the preprocessing stage, but you can launch it before the restoration process from the TOOLS > CROP ORIGINAL IMAGE menu. Once you have cropped your image during the guided restoration process you can not crop it again except after closing the image and reloading it again.

### Cropping an image in x y z

After you have verified your image parameters the Crop tool is launched if you press the YES button to the question 'Launch the cropper?'. The cropper will look as in Figure 13 (the image will be in gray scale mode if it is a single channel). Red lines indicate the borders of the proposed cropping region. This is computed from the image content and the microscopical parameters at launch time of the cropper.

The cropper allows manual adjustment of the proposed crop region. To adjust the crop region put the cursor *inside* the red boundary, press the left mouse button and keep it pressed to sweep out a volume. Accept the new borders by pressing the CROP button. Do not crop the object too tightly, because you would remove blur information relevant for deconvolution. Do not crop the image to make it too large along the optical axis Z, an aspect ratio close to 1:1:1 (or less than 1 for Z) is much better.

The three views shown are Maximum Intensity Projections (MIP's) along the main axes. The projections are computed by tracing parallel rays perpendicular to the projection plane through the data volume, each ray ending in a pixel of the projection image. The maximum intensity value found in each ray path is projected. For example, each pixel in the xy projection image corresponds with the maximum value in the vertical column of voxels above it.

By default the projections are over the whole dataset (including all the frames in time series), but this might be confusing sometimes. The small colored triangles can be used to restrict the projections within a specific range of slices.




Figure 12. Exporting a microscopical parameters template.

### Cropping an image in time

You also can reduce the number of time frames by selecting **TIME > SELECT FRAMES** from the **CROP** menu, as shown in Figure 14. This applies to time series (see *Time series* on page 8).

### Removing channels

You can remove channels from a multi-channel image (see p. 17) using the crop tool's channels selecting tool **CHANNELS > SELECT CHANNELS**.

### Stage 1: Parameter tuning

Stage 1 enables you to tune your parameter settings as set in the pre-processing stage. After you have finish cropping stage 1 is skipped and you will directly jump into the image histogram (stage 2). Still you may wish to tune your parameter settings afterwards. You can enter stage 1 by pressing the **RESTART** button in the latest stage.

### Stage 2: The image histogram

The next stage shows you the image histogram.

The histogram image is an important statistical tool for *inspecting* your image. It is included to let you spot problems that might have occurred during the recording. It has no image manipulation options as such, it just may prevent you from future recording problems.

The histogram shows the number of pixels as a function of the intensity (gray value) or groups of intensities. If your image is an 8-bit image (gray values from 0–255) the x-axis is the gray value and the y-axis is the number of pixels in the image with that gray value. If the image is more than 8-bits gray values are collected to form a 'bin' (for example gray values from 0–9 form bin 0, values from 10–19 form bin 1, etc.) The histogram is now the number of pixels in every bin.

The histogram from the demo image as shown in Figure 15 is of reasonable quality. The narrow peak you see at the left represents the background pixels, all with similar values. The height of the peak represents the amount of background pixels. Because in this particular image there are many voxels with a value in a narrow range around the background the peak is higher than the other.

In this case there is also a small black gap at the left of the histogram. This signifies an electronic offset (blacklevel) in the signal recording chain of the microscope (see page 40).

If a peak is visible at the extreme right hand side of the histogram it indicates saturation or 'clipping'. Clipping is caused by intensities above the maximum digital value available in your microscope. Usually, all values above the maximum value are replaced by the maximum value. On rare




Figure 13. The Crop tool.




Figure 14. Reducing the number of time frames.

occasions they are replaced by zeros. Clipping will have a negative effect on the results of deconvolution, especially with WF images. (See page 46).



The histogram stage is included for examining purpose only. It has no meaning for the deconvolution process that follows<sup>3</sup>.

### Stage 3: Estimate the average background in the image

In this stage the *average* background in a volume image is estimated. The average background is thought to correspond with the noise-free equivalent of the background in the measured (noisy) image. It is determined by searching the image first for a region with low values. Subsequently the value for the background is determined by searching in this region for the area with radius  $r$  which has the lowest average value. It is important for the search strategy that the microscopic parameters of the image are correct, in especially the sampling distance and the microscope type.

The following choices are possible here:

- **Lowest value** (Default): The image is searched for a 3D region with the lowest average value. The axial size of the region is around 0.3 micron; the lateral size is controlled by the *radius* parameter which is default set to 0.5 micron.
- **In/near object**: The neighborhood around the voxel with the highest value is searched for a planar region with the lowest average value. The size of the region is controlled by the *radius* parameter.
- **WF**: First the image is searched for a 3D region with the lowest values to ensure that the region with the least amount of blur contributions is found. Subsequently the background is determined by searching this region for the planar region with radius  $r$  that has the lowest value.

Press the ESTIMATE button to continue. You may now adapt the value if you like to, either by altering the value in the ESTIMATE BACKGROUND field or in the RELATIVE BACKGROUND field. Setting this last




Figure 15. **The image histogram.** The vertical mapping mode can be selected from linear, logarithmic or sigmoid.




Figure 16. **Estimating the average background.**

<sup>3</sup> Learn more about histograms in <http://support.svi.nl/wiki/ImageHistogram>



to -10, for example, lowers the estimated background with 10%. If you are done press ACCEPT to start the last stage, the true deconvolution process.

## Stage 4: Deconvolution

Before starting the actual iterative deconvolution run, stage 4 first carries out several pre-processing steps:

1. **Background:** This value is calculated in stage 3. You can verify whether this value represents areas in the image which you consider background by opening the Twin Slicer (see Figure 21) and moving the mouse pointer over the areas of interest. The voxel values are displayed below the image. Modify the value as you see fit.
2. A Point Spread Function (PSF, see below) was generated from the established microscopic parameters. This took place off the screen and is fully transparent to the user.
3. If the size of the computer's RAM is too small to deconvolve the image as a whole, it is split up in parts called 'bricks'. SVI's Fast Classic Maximum Likelihood Estimation (MLE) algorithm runs on the image or on all the bricks and fits the deconvolved bricks seamlessly together (see *Brick wise processing* on page 39).




Figure 17. Stage 4: starting the deconvolution.

## The Point Spread Function (PSF)

One of the basic concepts in image deconvolution is the Point Spread Function. The PSF of your microscope is the image which results from imaging a point object in the microscope. Due to wave diffraction a point object is imaged, 'spread out', into a fuzzy spot: the PSF. In fluorescence imaging the PSF completely determines the image formation. In other words: all microscopic imaging properties are packed into this 3D function. The PSF can be obtained in two different ways:

*1. Generating a theoretical PSF:* When a measured PSF is not available, Huygens Essential automatically uses a theoretical PSF. The PSF is computed from the microscopic parameters that come with your image and which you have double-checked in the pre-processing 'P' stage or in stage one. Because a theoretical PSF can be generated without any user intervention Huygens Essential does the calculation in the background without any notice.



Images affected by spherical aberration (S.A.) due to a refractive index mismatch are better restored through the use of theoretical depth-dependent PSF's. Read about S.A. on page 46, and how to correct it on page 48.

*2. Measuring a PSF:* By using the PSF Distiller or the tools in Huygens Professional you can derive a measured PSF from images of small (< 200 nm) fluorescent beads. You can load a previously measured PSF with FILE > OPEN PSF... (Main menu). If you load a PSF, Huygens Essential will automatically use it. If the measured PSF contains less

channels than the image, a theoretical PSF will be generated for the channels where there is no PSF available.

See *The PSF Distiller* on page 36 for more information on measuring a PSF.



A measured PSF should only be used for deconvolution if the image and the bead(s) were recorded with the same microscope at the same parameter settings as the bead image(s).

### The Classical Maximum Likelihood Estimation (CMLE) algorithm

Huygens Essential uses the Classical Maximum Likelihood Estimation for the deconvolution process. This method is an extremely versatile algorithm, applicable for all types of data sets<sup>4</sup>.

The following option values may be set:

#### Number of iterations

MLE-based deconvolution uses an iterative process that never stops if no stopping criterion is given. This stopping criterion can simply be a maximum number at which the process will stop. This value depends on the desired final quality of your image. For an initial run you can leave the value at its default. To achieve the best result you can increase this value.

Another stopping criterion is one based on the Quality change of the process, see *Quality threshold* below.

#### Signal to noise ratio

You have to make an estimation of the SNR from your recorded image. Inspect your image and decide if your image is noisy ( $\text{SNR} < 10$ ), has moderate noise ( $10 < \text{SNR} < 20$ ) or is a low-noise image ( $\text{SNR} > 20$ ). See *Signal to Noise Ratio (SNR)* on page 39, and some examples of noisy images in Figure 31.

#### Quality threshold

Beyond a certain amount of iterations, typically below 100, the difference between successive iterations becomes insignificant and progress grinds to a halt. Therefore it is a good idea to monitor progress with a quality measure, and to stop iterations when the change in quality drops below a threshold. At a high setting of this quality threshold (e.g. 0.1) the quality difference between subsequent iterations may drop below the threshold before the indicated maximum number of iterations has been completed. The smaller the threshold the larger the number of iterations which are completed and the higher the quality of restoration. Still, the extra quality gain becomes very small at higher iteration counts.

#### Iteration mode

In **FAST MODE** (highly recommended) the iteration steps are bigger than in **HIGH QUALITY mode**. More information can be read in the **DICTIONARY** from the **HELP** menu.

#### Bleaching correction

The data is inspected for bleaching. 3D images and time series of WF images will always be corrected. Confocal images can only be corrected if they are part of a time series, and when the bleaching over time shows exponential behavior.

<sup>4</sup> Huygens Professional also has Quick-MLE-time, Quick-Tikhonov-Miller, and Iterative Constrained Tikhonov-Miller.

## Stopping the MLE algorithm

Pressing **RESTORE** starts the iterative MLE algorithm; a **STOP** button appears. Pressing **STOP** halts the iterations and retrieves the result from the previous iteration. If the first iteration is not yet complete a empty image will result.

## Finishing or restarting a deconvolution run

When a deconvolution run is finished use the *Twin Slicer* (page 21) to inspect the result in detail. Depending on the outcome of that you can choose **AGAIN**, **RESUME** or **ACCEPT**:

**AGAIN** discards the present result, and re-runs the deconvolution, possibly with different parameters.

**RESUME** re-runs the MLE procedure without discarding the result, and with the possibility to change the deconvolution parameters. The software will ask you to continue were you left off (keeping improving the image, quite recommendable) or to start from the raw image again. A new result will be generated to compare with the previous one, for instance using the *Twin Slicer*. You can repeat this several times.

**ACCEPT** proceeds to the final stage or, if the data was multi-channel, to the next channel (see page 17). If you generated several results by resuming the deconvolution you will be asked to select the best result as the final one, that will be renamed to 'deconvolved'. The other results will remain as well in case you want to save them.

## z-drift corrector for time series

For 3D time series the program pops-up an additional tool that enables you to correct for movement in the *z* (axial) direction that could have been occurred for instance by thermal drift of the microscope table. In case of a multi channel image (see p. 17), the corrector can survey **ALL CHANNELS** and determine the mean *z* position of the channels, or it can take **ONE CHANNEL** as set by the **REFERENCE CHANNEL** parameter.

After determining the *z* positions per frame, the *z*-positions can be filtered with a **MEDIAN**, **GAUSSIAN** or **KUWAHARA** filter of variable width. When the drift is gradual, a gaussian filter is probably best. In case of a drift with sudden reversals or outliers a median filter is best. In case the *z* positions show sudden jumps, we recommend the edge preserving Kuwahara filter.




Figure 18. The *z*-drift corrector.

## Saving the result

### Saving the restored image

After each deconvolution run you can save the result. Select the image to be saved and do **FILE > SAVE 'IMAGE NAME' AS...** in the menu bar. You can save the image as an Image Cytometry Standard (ICS or ICS2) image file, a TIFF file series, an Imaris-Classic file or a Biorad .pic file. Only the ICS and ICS2 file type save all the microscopic parameters. For information on how to proceed with multi-channel data see *Joining the results* on page 17.



The ICS file format actually uses two separate files: a header file with `.ics` extension and other, much bigger and with the actual image data, with `.ids` extension. On the other hand the newer ICS2 file format uses only one single `.ics` file with both the header and the data together.


### Saving your Task report

Select from the Main menu bar `TASK > SAVE TASK REPORT` to store the information as displayed in your *Task report* tab deck.

## Multi channel images

Fluorescent microscopes can usually register different emission wavelengths (almost) simultaneously, allowing you to image different dyes on the sample. In the terminology of the Huygens Software, one channel in an image refers to the intensity distribution recorded at a given fixed wavelength, independently of what device made the acquisition. Thus, it is a logical channel of stored data, and not necessarily a physical channel (as all the image channels could have been measured by a single photomultiplier, for instance).

Multi channel images can be deconvolved in a semi automatic fashion, giving you the opportunity to fine tune the results obtained with each individual channel. *After the preprocessing stage* the multi channel image is split into single channel images named 'channel-0', 'channel-1', and so on. The first of these is automatically selected for deconvolution. To deconvolve it, proceed as follows:



**Deconvolving a channel in a multi channel image** Figure 19. *Deconvolving a two channel image.*

The procedure to deconvolve a channel in a multi channel data set is exactly the same as for a single channel image. You can therefore do multiple reruns on the channel at hand, just as you can with single channel data. When you are done press `ACCEPT` in the last (stage 4) screen. This will cause the next channel to be selected for restoration. Proceed as usual with that channel and the remaining channels. If you do not want to process all the channels in an image you may skip one or more channels.

### Joining the results

When you press `ACCEPT` for the last channel you enter a screen which allows you to select the results which you want to combine into the final deconvolved multi channel image. This means that up to this point you can still change your mind as to which of the results you want to combine, even in what order. Once you press `ACCEPT` a multi channel image named *Restored* is created. To save it go to the `FILE > SAVE 'RESTORED'` AS menu.

## Batch processor

Once you know how to deal with a particular kind of dataset and are sure of the restoration parameters (see *Improve the estimated parameters* on page 48), you can restore a couple or more of similar datasets automatically. This is called **batch processing**.

A batch process is made of independent image restoration tasks (one per image) that are executed one by one until all are finished. Depending on the multithread capabilities of the computer multiple tasks can be executed in parallel.

You can for example program batch scripts using Huygens Scripting, that enables you to run scripts written in Tcl, using the extensive set of Tcl-Huygens image processing commands.

You can also configure batch processes easily using the interactive Huygens Batch Processor. The Batch Processor is the tool to do large scale deconvolution of multiple images within the Huygens Essential.

### Start up

To launch the Batch Processor first open Huygens Essential, then click on the menu DECONVOLUTION > BATCH PROCESSOR.

### Window description

These are the different elements that form the window:




Figure 20: Batch Processor main window.

- **Save location** is the location where the resulting images will be placed during the batch run. With the two folder buttons you can respectively select a location or create a new one in the currently selected folder.

- **Current task lines** shows a list of tasks (empty at start). Tasks are jobs that will be processed by the Batch Processor one by one. Each task line consists of an image, a microscopic template and a deconvolution template. These templates can be updated after a task line is added to the list to tune the values in each particular case.
- **Add Task Line(s)** to add tasks to the list, one by one or many at the same time.
- The text tabs (**Info**, **Progress**, and **Report**) give information about the whole process in its different stages.
- **Run** button, to launch the process when everything is configured.

## Usage

### Adding one task

First let us try adding a single task:



- Click on **SELECT AN IMAGE FILE...** and browse to select one raw image that you want to deconvolve. End by clicking **OPEN**.
- Then click on **MICRO TEMPLATE...** to load on the microscopic parameter set that you use in Huygens Essential. If none exists, you must **CREATE** one. The **CREATE** actually opens a template editor, where you can load and save many different templates (overwriting the existing ones, or creating new ones).
  - Tip: prepare your templates in the Huygens Essential main window before selecting them in the Batch Processor (see page 10). This will allow you to use the parameters saved inside an image.
- Repeat the same procedure with the **DECON TEMPLATE...** to establish a set of restoration parameters. In the Batch Processor, apart from the wizard CMLE restoration algorithm, you can also use the QMLE, a quicker version that works very well with low noise images<sup>5</sup>.
- Finally, the icon with the single green plus sign will be enabled: click on it to add the task.

You will see the task included in the list of tasks. By now it is the only one, we will add more later. You can now click on the different task elements:

- the image name, to change the data file the task is executed with.
- the microscopic template, to change the microscopic parameters You can modify them or load a new template.
- the deconvolution template, to change the restoration parameters. You can modify them or load a new template.
- the "trash can" icon, to remove the task from the list.

### Adding multiple tasks

Now that you have learnt how to add one task to the job list and explored its capabilities, we can speed things up and add many tasks at the same time.

<sup>5</sup> See <http://support.svi.nl/wiki/RestorationMethod> for more information.

You can use the button with the two plus signs for that:



That will open a new dialog, where you can select multiple images that will be deconvolved with the same microscopic and restoration templates, as specified.

You can add all image files present in a **directory**, or add multiple **files**. In the second alternative, multiple files can be selected by holding down the control key (apple in the Mac). Once you have accepted a list with all the images, select the templates to be used with them, and click **ADD** to generate all the task lines. All of them will use the selected templates.

Repeat the procedure to add more image or use different templates. Click **CANCEL** to finish.

When you change one task's template you see that the template name is also changed, indicating that some parameters were modified locally. These changes do not affect the template as it is used by the other tasks.

#### **Running the batch job**

Make sure you selected the desired location to save the results.

When you have your batch process configured, you can save it for future reference: in the menu, do **FILE > SAVE**.

When everything is ready you just click the **RUN** button to start the jobs.

The progress of the Batch Processor and the report for each individual task are shown on the text tabs of the *Batch Processing Overview*. The status of each task on the task list changes accordingly to the evolution of the process.

The restored images are saved on the selected destination directory as soon as they are ready.

## **Menus**

### **File menu**

By using this menu you can **CLOSE** the current process, create a **NEW** task list, **OPEN** previously saved processes and **SAVE** them to the disk, or **APPEND** other task lists to the current one.

You can also **EXIT** the Processor.

### **Task menu**

By using this menu you can **RUN** the tasks, and **SAVE BATCH PROGRESS** or **SAVE BATCH REPORT** to text files on the disk.

## **More information**

Updated information about the Batch Processor can be found in <http://support.svi.nl/wiki/BatchProcessor>.

## 4. Huygens Essential Visualization Tools


Huygens Essential provides different tools for data visualization.

### The Twin Slicer

This allows you to compare a deconvolution result with the original, but also different deconvolution results obtained from the same original.

Open the Twin Slicer by double clicking on one of the thumbnail images in the main window (or by right-clicking on the image thumbnail, then **SHOW ON TWIN SLICER**). The Twin Slicer will show the selected image on the left. By clicking on the menu bar below the image you can select a different image. Likewise, the bar below the right display field gives access to one of the other images currently present in Huygens Essential (see Figure 21). Currently only two images with the same dimensions can be displayed at the same time.

If you select the same image twice, you may compare different slices from the same image. For this, first move the slider until the desired position, then click on one of the image's name to **DISABLE** the action of the slider on it: the slider will then affect only the other image view (see Figure 22).



Click the name button to select an image, disable the slider function or center the image.

Figure 21. **The Twin Slicer**, applied to a 3D dataset. The selected display settings are highlighted in blue.

Pixel intensity values for the cursor position on the image are displayed at the bottom of the window. You can move the image by clicking the left-mouse button and keeping it pressed while moving the image to the desired position. You can center the image again by selecting **CENTER IMAGE** from the name button.

With the slider you can slice your images along the three axis, depending on what is the selected slice view. You can also change the zoom factor, the color mode and the contrast mapping mode.

### Color mode

- Gray: pixel values are assigned different shades of gray ranging from black for the lowest values to white for the highest values.

- False color: pixels values are assigned different colors ranging from black/dark purple for the lowest values to bright red for the highest value.
- True color: pixel values are assigned different shades of a particular color ranging from black for the lowest values to the brightest possible shade for the highest values.




Figure 22. The Twin Slicer showing two slices of the same 3D image.

Multi channel images are always rendered with a true color scheme, otherwise the information from the different channels will result in very confusing images.

### Contrast mapping mode

- Linear (default): In this mode the pixel values are mapped to screen buffer color intensities in a linear fashion. Note that the actual translation of the screen buffer values to the actual brightness of a screen pixel is usually quite non-linear.
- Compress: Where an image contains a few very bright spots and some larger darker structures using Linear mode will result in poor visibility of the darker structures. Restoration of such images is likely to further increase the dynamic range resulting in the large structures becoming even dimmer. In such cases use the compress display mode to increase the contrast of the low valued regions and reduce the contrast of the high-valued regions. Another way to improve the visibility of dark structures is the usage of false colors (see above).
- Widefield (WF) mode: In restoring widefield images it sometimes happens that blur removal is not perfect, for instance when one is forced to use a theoretical point spread function in sub optimal optical conditions. In such cases the visibility of blur remnants can be effectively suppressed.

### Time series

If you open the Twin Slicer on a time series a second slider is added. Both the time-slider and the spatial-slider have a swing option. When the spatial swing is pressed the slider moves back and forth; when the time swing button is pressed the slider only moves forward, i.e. in the positive time direction.




Figure 23. The Twin Slicer with time slider.



## The MIP Renderer

The Maximum Intensity Projection (MIP) Renderer is part of the Huygens Essential since release 2.6.0p4 and enables you to obtain a spatial projection of your data from a given viewpoint (see Figure 24).

The renderer projects, in the visualization plane, the voxels with maximum intensity that fall in the way of parallel rays traced from the viewpoint to the plane of projection. Notice that this implies that two MIP renderings from opposite viewpoints show symmetrical images.

To start the MIP Renderer, right-click on an image's thumbnail to open the contextual menu, then select **SHOW IN MIP RENDERER**.

Select your viewpoint by moving the 'Tilt' and 'Twist' sliders, or by dragging the mouse pointer on the large view (that will be empty at first, before your first rendering). Also try changing the zoom. You will see how the preview thumbnail changes. When you have set all the rendering options, click **Render** to create the final view, that you can save as a TIFF image.




Figure 24. *The Maximum Intensity Projection (MIP) Renderer*

You will find different rendering options on the window, and also in the Options menu. The configurable parameters are the rendering size and quality, the appearance of the bounding box, and the mode of the soft thresholds applied to the image channels.

### Soft threshold

A soft threshold is a preprocessing tool that reduces the background in the image, so voxels with intensity values below the threshold value become more transparent. Contrary to a standard threshold, that is 'all or nothing' (values above the threshold are kept, values below it are deleted), the soft threshold function handles images in a different way. It makes smooth transitions between the original and the deleted values. If the original value  $S$  of a voxel is  $S > (\text{threshold value} + \text{range}/2)$  then the final filtered value  $D$  does not change ( $D = S$ ). If  $S < (\text{threshold value} - \text{range}/2)$  then the voxel is 'deleted' ( $D = 0$ ). For the values in between, a smooth function is applied: if  $(\text{threshold value} - \text{range}/2) < S < (\text{threshold value} + \text{range}/2)$  then  $D = f(S)$  according to a shape function, which in this case is a sinusoidal. By changing the parameters in **OPTIONS > SOFT THRESHOLD** mode from **HARD** to **SOFT**, you progressively increase the 'range' value, thus broadening the transition from the original to the deleted values.

You can apply different soft thresholds to the different image channels.

### Rendering a movie

With the MIP Renderer you can also make an animation of your image, changing the viewpoint in different frames. Select the viewpoint coordinates for the first frame, then click **SET > HOME**. Select now the viewpoint coordinates for the last frame, and click **SET > END**. (You can now go to the last or the first frame by clicking **Go > END** or **Go > HOME**). Select all the rendering parameters, including the total number of rendered frames for the movie (**OPTIONS > ANIMATION FRAME COUNT**). Finally, click **ANIMATE**, and select a directory to save the TIFF frames to. You can later load and edit

these TIFF images with your favorite animation tool. For instance, you can use the `convert` tool from ImageMagick (<http://www.imagemagick.org>) to make a GIF animation, using

```
convert -delay 20 animatedMip*.tiff animatedMip.gif
```

You can now place this single file GIF animation directly on your web page, as most of the Internet browsers currently available can handle this kind of movie files.

With the appropriate codec, you can also use `convert` to make a MPEG animation. See the ImageMagick website for more details.

If your image is a time series, you can also make an animation along time frames.

## The SFP Renderer

Starting from Essential 2.5, a simplified version of SVI's high end volume renderer FluVR (Fluorescence Volume Renderer) is available for visualizing your volumetric object from a selectable viewpoint. Like FluVR, this renderer is based on taking the volume image as a distribution of fluorescent material, simulating what happens if the material is excited and how the subsequently emitted light travels to the observer. The computational work is done by the Simulated Fluorescence Process (SFP) algorithm.

The ray-tracing technique does not require a special graphical board as the polygon based techniques do.

To start the SFP renderer, right-click on an image's thumbnail to open the contextual menu, then select `SHOW IN SFP RENDERER` (see Figure 25).

### Summary

A virtual light source produces excitation light that illuminates the object. This casts shadows either on parts of the object itself or on a table below it. The interaction between the excitation light and the object provokes the emission light, that also interacts with the object before it reaches the eye of the viewer (see Figure 26).

### SFP fundamentals

The voxel values in the image are taken as the density of a fluorescent material. If the voxels are multiparameter ("multi-channel" in microscopic parlance, see page 17) each parameter is taken as a different fluorescent dye. Each dye has its specific excitation and emission wavelength with corresponding distinct absorption properties. The absorption properties can be controlled by the user. The different emission wavelengths give each dye its specific color.





Figure 25. **The SFP volume renderer.** Up right the preview subsampled image that acts directly on the slider position. The actual rendering is started by pressing the **RENDER** button.



To excite the fluorescent matter light must traverse other matter. The resulting attenuation of the excitation light will cause objects, which are hidden from the light source by other objects, to be weakly illuminated, if at all. The attenuation of the excitation light will be visible as shadows on other objects. To optimally use the depth perception cues generated by these shadows a homogeneous plane (the gray table) below the data volume is placed on which the cast shadows become clearly visible.

After excitation the fluorescent matter will emit light at a longer wavelength. Since this emitted light has changed wavelength it is not capable to re-excite the same fluorescent matter: multiple scattering does not occur. Thus only the light emitted in the direction of the viewer, either directly or by way of the semi reflecting table is of importance. By simulating the propagation of the emitted light through the matter, the algorithm computes the final intensities of all wavelengths (the spectrum) of the light reaching the viewpoint. By default the first channel (ch-0) is the red object, the second channel (ch-1) is the green object and the third (ch-2) is the blue object.



The properties of the interaction between object and light (transparency), both for excitation and emission, as well as the viewpoint, can be adapted interactively by the user to produce different sceneries. Since the volume rendering process is rather compute intensive, a preview image is displayed (see Figure 25). Apart from the viewpoint settings and the optional zooming, the following sliders affect the image:

- Transparencies:
  - Excitation: The transparency of the object for the excitation light. The less transparency the more shadow is casted on the subsequent voxels and on the table.
  - Emission: The transparency of the object for the emission light. The lower the transparency for the emission light the more difficult it is to peer inside or trough the object.
- The Characteristic object size affects both the excitation and the emission transparency. While traveling through the object, the light intensity is attenuated to some degree. This enables us to define some definition for penetration depth at which the light intensity is decreased to some extent, say 10% of its initial value. This penetration depth should be in line with the object size. A transparent object is small with respect to the penetration depth. Thus for the same physical properties of the light one object can be transparent while the other is oblique due to its size. To find a reasonable range in transparencies the object size may be altered. At start-up the object size is computed from the microscopic sampling sizes and number of pixels the image is composed off. If your image has not the

**Figure 26. With the SFP renderer excitation and subsequent emission of light of fluorescent materials is simulated.** Each subsequent voxel in the light beam (excitation) is affected by shadowing from its predecessors. The transparency of the object for the emission light controls to what extent the viewer can peer inside the object.

The light source is drawn here inside the figure, but in real is placed at infinite distance as to make the light rays parallel.

The renderer in Essential is in non-perspective mode (so called orthogonal projection), i.e. the viewpoint is at infinity.

correct parameters (for example a TIFF series) the object size is set according to the default parameters as set by the Huygens Essential software and may not be related to the actual object size.

- Frame: Time series (frames) can be handled. For a 3-D image this slider is inactive and set to Frame 0.
- Object brightness: Intensity of the virtual light source.
- Soft threshold: Preprocessing tool that reduces the background in the image. Voxels with gray values below the threshold value become more transparent. See a detailed explanation in the MIP Renderer section, on page 23.

Use the **RENDER** button to start the actual rendering. The result can be saved as a TIFF-image (**FILE > SAVE**).

Use **LIGHT DIRECTION** to alter the direction of the light.

You can open as many SFP windows as you like.

### Rendering a movie

With the SFP Renderer you can also make an animation of your image, changing the viewpoint or the time coordinate in different frames. The procedure is analogous to the one explained for the MIP renderer on page 23.

---

## The Surface Renderer

The Surface Renderer is available from the Huygens Essential version 3.0 onwards, and enables you to represent your data in a convenient way to clearly see separated volumes. Because this Surface Renderer is based on fast raytracers, there is no need for any special graphic card as would be necessary for conventional polygon based techniques.



The Surface renderer is an extended optional tool, and is enabled by a *v* flag in the license string (see *License string details* on page 49).

To start the Surface Renderer, right-click on an image's thumbnail to open the contextual menu, then select **SHOW IN SURFACE RENDERER**. Let the renderer initialize.

You can find three graphic pipes to redirect your image data channels to: two surface pipes, and one MIP pipe. These can be activated independently.

Use the threshold slider to apply different thresholds to your data channels, to select what voxels are considered to shape volumes. Connected voxels after the threshold determine independent volumes, that will be represented by the 3D isosurfaces containing them, with different colors. You can use up to two different data channels for surface rendering, one in each of the two available surface graphic pipes. The color of the different objects inside a channel can be modified with a selector (see *Hue Selector* below).

There is a third graphic pipe to redirect data to the rendered image: the MIP pipe works projecting the voxels with maximum intensity laying in the path of the rays traced along the viewing direction (see *The MIP Renderer* on page 23). In combination with the surface pipes, you can obtain very clear representations of the different objects in your image. The MIP rendering of one channel can be a good spatial reference for the objects in the other channels.





Figure 27. *The Surface Renderer.*

You can control the transparency and the brightness of the rendered surfaces with the corresponding sliders, independently in each graphic pipe.

Select the viewpoint by moving the 'Tilt' and 'Twist' sliders, or by dragging the mouse pointer on the large view. Also try changing the zoom. When you have your rendering ready, save it to a TIFF file with **FILE > SAVE**.

Other available options to change are render size and transparency depth, accessible through the **OPTIONS** menu.

### Hue Selector



The hue selector is a component that allows you to select the color range (actually, the hue property of the color) in which the different objects of each channel are displayed by the Surface Renderer or the Colocalization Analyzer. Thus, objects belonging to different channels can be represented with very different hue ranges to make them clearly distinct, but also with some gradual differences inside the selected range to distinguish independent objects. You can also collapse a range to have all objects in a channel displayed with exactly the same color.

### Transparency depth

This option controls how different surfaces are seen through the others. The effects are mainly visible when you have some objects intersecting with other. With the **SIMPLE** depth, only the piece of surface closest to the viewer's eye screens the others behind it, with its corresponding transparency level.

With the **NORMAL** depth, up to two pieces of surface are considered to screen other objects. Thus, one object *B* *inside* the surface *A* will appear less screened than a third object *C* *behind* *A*. *B* is only screened by the piece of surface *A* closer to the viewer's eye, while the object *C* is screened by two pieces of surface *A*.

The **DEEP** option will consider many more screening levels, making the final rendering more complex.

---

## 5. Huygens Essential Analysis Tools

Huygens Essential is extended with new tools for interactive analysis of 3D and 4D microscopic images. The Colocalization Analyzer is present in Huygens Essential since version 3.0, and the Object Analyzer since version 3.1.

---

### The Object Analyzer

The interactive Object Analyzer tool allows you to obtain statistics of individual objects by clicking on them, or analyzing all objects with a single button press.

In this context, an **object** is a distinct group of interesting voxels that are spatially connected one to another. Interesting voxels are distinguished from the background by using a seed and threshold criterion.

Therefore defining an image's objects implies

1. **segmentation**: separating interesting voxels from the background according to a given criterion, and
2. **labeling**: grouping them together.

This is done interactively by the Object Analyzer. To remove too small objects in an early stage from the analysis a garbage level can be set below which objects are discarded. After that, detected objects are automatically labeled and sent to a continuous isosurface renderer.

An **isosurface** is a 3D surface representation of points with equal values in a 3D data distribution. Is the 3D equivalent of a 2D contour line.



The Object Analyzer is an extended optional tool, and is enabled by a A flag in the license string (see *License string details* on page 49).

### How to use the Object Analyzer

To start the Object Analyzer in Huygens Essential, right-click on an image's thumbnail to open the contextual menu, then select **SHOW IN OBJECT ANALYZER**. Alternatively, you can select an image's thumbnail, then in the menu bar select **VISUALIZATION > LAUNCH OBJECT ANALYZER**. Wait for the analyzer to initialize and to compute the first image based on the default parameters. It is recommended to use a deconvolved image.




Figure 28. *The Object Analyzer*, with segmented objects represented as isosurfaces. An area is selected to be extruded and define a 3D ROI.

### Mouse mode

Depending on the selected mode, the mouse has different behaviors:

- **Rotate object** interacts with the full image to rotate it in the space, by dragging the mouse pointer on the object view. That can also be achieved by moving the 'Tilt' and 'Twist' sliders.
- **Measure object** lets you click on different defined objects and obtain the local statistics.
- **Define ROI** lets you define different regions of interest (ROI), discarding objects out of them.
- **Delete object** lets you delete irrelevant objects one by one. Both this and the ROI selector act only on the current window, the original image is not changed.
- **Set measurement anchor** lets you select and unselect objects to be 'anchors', references for distances to other objects when asking for local statistics. It is possible to select a group of anchor objects. The location of the anchor is then the centre of mass of all anchor objects.

Independently of the mouse mode, you can at any time use the right mouse button (middle in the Mac) to drag the image and **recenter** the objects.

The different window elements are explained with a short text tip on the bottom of the window when the mouse passes over them.

First we have to select what is the data to be analyzed.

### Frame selector

In time series, you can select the interesting time frame with this slider.

## Render pipes

You can find three graphic pipes to redirect your image data channels to: two isosurface pipes, and one MIP pipe. These can be activated independently, by selecting a channel data from the image. Objects, as defined above, are represented in the isosurface pipes. The MIP pipe is used as a reference.

You can use up to two different data channels for isosurface rendering, one in each of the two available surface graphic pipes.

Use the threshold sliders to apply different threshold levels to your pipes, to select what voxels are considered to form objects based on their intensities. The remaining adjacent voxels are grouped together to determine independent volumes, that will be represented with different colors:

1. The image is segmented into objects by an effective **seed and threshold** level and 6-connection technique. The image is first explored for features that have intensities above the seed level, then the selection is expanded around this seed to capture all the neighbor voxels down to the threshold. If the seed is set to zero, then the segmentation is done by a standard threshold level<sup>6</sup>. Voxels connected by one face are grouped together.
2. To remove too small objects in an early stage from the analysis a **garbage volume** (in pixels) can be set below which objects are discarded.
3. The segmented and automatically labeled image is shown as a colored isosurface image. Each object is assigned a different hue from a range of hues selectable with the hue selector (see p. 27).

Because the isosurface renderer can recompute a new isosurface rapidly, the thresholds can be manipulated interactively.

## Object removal

Uninteresting objects can be deleted one by one by clicking on them in the *delete object* mouse mode explained above. All the voxels in a channel inside the removed object are not considered anymore, even when you change the thresholds again.

## Regions of interest (ROI)

You can define a region of interest (ROI) in the data, to discard objects out of it for further analysis.

1. Select the *define ROI* mouse mode.
2. Draw an area on the isosurface pane keeping the mouse button pressed down. When you release it, the area selection will be automatically closed. The selected area is extruded perpendicular to the screen to create a 3D selection.
3. Select on what data channel you want to operate. Your options are
  1. One channel, as selected in the Primary graphic pipe.
  2. All channels.
4. Click on the appropriate button what do you want to delete: objects **INSIDE** or **OUTSIDE** the selected area.

You will see all the excluded objects disappear and the volume enclosed in the ROI become highlighted. The way the ROI is represented is controlled by another graphics pipe, whose controls are grouped in the *ROI surface* pane. As you can define different ROIs for different channels, you can use different visualization parameters for each of them. Switch the highlighting off to leave just the data.

---

<sup>6</sup> For more explanations see <http://support.svi.nl/wiki/SeedAndThreshold>

You can repeat the select-delete procedure in this or another view to refine your ROI definition as much as you want.

To recover the original data again, press the **RELOAD** button.

The **EXPORT** button takes the image data inside the currently defined ROI and exports it to overwrite the original image in Huygens Essential, that was not modified up to this moment. Data outside the ROI is removed from the image, that can now be saved in a new file. This makes the Object Analyzer a powerful cropping tool, useful as a independent preprocessing tool for visualization or for example colocalization analysis.

The **RELOAD** button always reloads the image in the main window. That means that it will reload the result of your last Export, if you did any.

### Data analysis

When the *measure object* mouse mode is active, local statistics are returned in the lower window pane when clicking on different objects. As these properties (and the number of objects) vary with the selected thresholds, first a report of the current conditions is shown:

```
Conditions: pipe 0 chan 1: thresh 136.7 seed 170.8 garb. 1 objects 4
           pipe 1 chan 0: thresh 256.1 seed 320.1 garb. 1 objects 4
```

That reports the active isosurface pipes, what data channels they are using, the current seed and threshold levels (in absolute values), the garbage level in pixels, and the total number of found objects per pipe.

**Notice** that the object statistics pane is a text box. Output will be placed in the current cursor position. Make sure to put it at the end again after selecting text for copying.

The following statistics are reported for each object:

- Object reference: Surface pipe (primary or secondary), position (front or back), data channel and label number
- Voxel count
- Volume enclosed by the local isosurface
- Surface of the local isosurface
- Sum, minimum and maximum of all object voxels intensities
- Length of the object
- Caliper width in the axial and lateral direction of the microscope
- Center of mass of the object

It may happen that clicking at one point of the screen produces an ambiguous selection of objects, because to have two of them, one behind the other, belonging to different graphic pipes. That will make the analyzer to report information about the two pointed objects. In cases like that is when the position information (front or back) is more relevant. There is no ambiguity when the two overlapping objects belong to the same pipe: in that case information is always reported about the object closer to the viewer. If they belong to the same pipe the overlapping is only apparent and due to the perspective projection: distinct objects inside a pipe are never in contact in space.

There is an **ANALYZE ALL** button to obtain information about all the presently defined objects on the primary pipe, without having to click on them one by one.

The width of the objects is measured using a virtual caliper held perpendicular to the length axis. The caliper is moved along the length axis to find the largest width while being rotated around.



Because microscopic data, even deconvolved microscopic data, often shows orientation dependent imaging due to the lower axial resolution, a distinction is made between width measurements where the measurement direction is in the axial direction vs. the lateral direction. The caliper rotates 180 ° in space along the length axis to find maximum widths, but this rotation is split in two 90 ° ranges: one where the component along the optical axis is dominant, and another one where lateral components are larger. Both values are reported.

When an **anchor** object was previously defined (in the *set measurement anchor* mouse mode), distance to the anchor is also reported at the top of the object statistics pane for the last pointed object.

It is possible to select a group of anchor objects. The location of the anchor is then the center of mass of all anchor objects. A button is available at the statistics pane to remove all the anchor positions.

### Visualization parameters

You can select the **viewpoint** of the rendering by moving the **TILT** and **TWIST** sliders, or by dragging the mouse pointer on the object view when the *rotate object* mouse mode is active. Also try changing the **Zoom** slider.

You can control the **transparency** and the **brightness** of the isosurface pipes independently with the correspondent sliders, or just set them off.

The **color range** in which the objects are displayed can be modified using the hue selectors (see p. 27).

Other available options to change are **render size** and **transparency depth** (see p. 27), accessible through the Options menu. The render size can be set to fit the size of the canvas, which is the pane where the rendering is displayed (and whose size depends on the size of the window)

### Save the results

All the obtained information can be saved to external data files through the **FILE** menu.

The exported image, as viewed on the screen, will have the size determined by the *render size* visualization parameter mentioned above. Notice that if the render size is too big it might not fit the display pane, but it will still be saved as rendered. You can shift the rendering by dragging it with the middle mouse button.

The exported text file includes the calculated object statistics as shown in the lower data pane.

### Read more

More information is available at <http://support.svi.nl/wiki/ObjectAnalyzer>.

---

---

## The Colocalization Analyzer

The Colocalization Analyzer provides information about the amount of spatial overlap between structures in different data channels, for 3D images and time series.

As this overlapping can be defined in many ways, Huygens gives you the colocalization coefficients most commonly used in literature: Pearson, Overlap, and Manders  $M$  and  $K$  (see *Colocalization Theory* in the SVI-wiki<sup>7</sup>).



The Colocalization Analyzer is an extended optional tool, and is enabled by a C flag in the license string (see *License string details* on page 49).

## Iso-colocalization object analysis

One of the features of the colocalization analyzer is iso-colocalization object analysis. It allows you to quickly determine the properties of the different colocalization regions in your data. This is realized by visualizing the colocalization map as iso-colocalization surfaces. In this way regions in which the degree of colocalization exceeds a certain value become objects. By clicking on the objects local colocalization parameters are computed and reported. To relate the iso-colocalization objects to the original data the surface objects can be blended with a MIP projection of the data (p. 23). The color range in which these objects are displayed can be modified using a hue selector (p. 27).

## How to use the Colocalization Analyzer

To **start** the Colocalization Analyzer in Huygens Essential, right-click on an image's thumbnail to open the contextual menu, then select **SHOW IN COLOCALIZATION ANALYZER**. Alternatively, you can select an image's thumbnail, then in the menu bar select **VISUALIZATION > LAUNCH COLOCALIZATION ANALYZER**. The image must be multi channel (see p. 17) as the colocalization is based on the overlapping of different channel intensities. Wait for the analyzer to initialize and to compute the first MIP projection of the data.

You can select the viewpoint of the MIP projection by moving the **TILT** and **TWIST** slides, or by dragging the mouse pointer on the object view. Also try changing the zoom.

At the beginning the iso-colocalization surface sliders are deactivated, as we have to calculate a colocalization map first.

First we **select the data** to analyze. For time series, the "frame selector" selects the time coordinate.

We follow the usual naming convention in colocalization theory for the two compared channels: Red (R) for the first channel, Green (G) for the second channel. We can select, in the lower part of the window, which data channels from our image are the Red and Green channels to be compared.

A two-channel **histogram** is calculated by default, and updated whenever we change the Red or Green channel selection. This histogram is already an indication of the degree of overlapping between the selected channels: for two channels with a high degree of overlapping, the histogram pixels trend to concentrate along the diagonal  $y = x$  line. On the other side, total absence of overlapping would produce a 2D histogram with values only on the coordinate axes<sup>8</sup>.

By moving the colored **background** lines in the histogram, or by changing the numeric values in the input fields, we specify the backgrounds for the two selected channels. These values are subtracted from the voxels intensities when calculating the coefficients (if the result is negative, it is understood as a zero). Generally the colocalization coefficients depend much on correct estimation of the image background and resolution. For these reasons we strongly recommend to compute colocalization coefficients only on deconvolved images.

<sup>7</sup> <http://support.svi.nl/wiki/ColocalizationTheory>

<sup>8</sup> Read more on <http://support.svi.nl/wiki/TwoChannelHistogram>




Figure 29. *Colocalization Analyzer* with surfaces for regions of 'overlapping'.

Backgrounds are for removing signal that would lead to spurious colocalization, prior to the calculation. This is intended for minor tuning, or just in case you really need to calculate colocalization in raw images that still have all the measured background. See *Backgrounds vs. thresholds in colocalization* below.

Next we select what **colocalization map** we want to calculate: Pearson, Overlap, Manders M, K, or H.

Notice the difference between **maps and coefficients**: the colocalization coefficients parametrize the degree of colocalization of the full image, while a colocalization map parametrizes the colocalization locally. In a map, a single colocalization value is calculated per voxel creating a 3D distribution that is represented in a 3D image by isosurfaces.

The Colocalization Analyzer computes only the map selected by the user, but it always computes all the available coefficients.

Finally we click on **COMPUTE**.

The *Frame colocalization coefficients* pane will show all available colocalization coefficients for the whole frame.

The obtained colocalization map is represented in the renderer window by iso-colocalization surfaces: regions in which the degree of colocalization exceeds a certain value become objects. This

"certain value" can be controlled by the threshold slider in the iso-colocalization surface parameters. You can control the transparency and the brightness of this surface pipe with the correspondent sliders, or just set this pipe off. The color range in which the objects are displayed can be modified using a hue selector (see page 27).

Some modes generate two-channel colocalization maps: colocalization of Red respect to Green, and vice versa, e.g. in case of the Manders M1 and M2 coefficients. In these cases, the iso-colocalization surface parameters will offer the possibility of rendering any of the two channels, and thus the threshold is referred to the active one.

By clicking on the rendered objects local colocalization parameters are computed and reported.

There is one Maximum Intensity Projection (MIP) pipe available to redirect your data channels to. The MIP rendering of one original channel (maybe none of the ones used for colocalization) or the two channels under analysis can be a good spatial reference for the objects from the colocalization map. When you select an original channel, you can use the threshold slider to select what voxels are considered for the MIP rendering, depending on their intensities. Notice that here the threshold is simply used for representation. If you select both R and G channels to be rendered, their correspondent backgrounds as selected in the histogram will be used as projection thresholds. As with the surface pipes you can also control the transparency and the brightness of this MIP.

All the obtained information can be saved to external data (text or image) files through the FILE menu.

### Backgrounds vs. thresholds in colocalization

**Backgrounds** are for removing signal prior to the calculation. In an ideally restored image that would not be necessary, because all the signal present in the image is "good signal". The background would have been removed during the image restoration.

**Thresholds** are to split the colocalization maps in two regions: what are interesting objects, and what are not. Local colocalization values are calculated for every image voxel, but the zero value would be very rarely achieved. Most of the times you have some non-zero colocalization everywhere, but you are probably interested in regions where colocalization exceeds, let us say, a value of 20 % of the maximum.

With the backgrounds, you discard data based on voxels intensities: high intensity regions are left, low intensity regions are disregarded. But that says nothing about colocalization! You can very well have high colocalization levels in regions with low voxel intensity.

Therefore you should not use backgrounds to remove signal that can still have some colocalization level. Ideally you consider all the signal, then study the colocalization levels using thresholds to split regions of high colocalization from regions with low colocalization, something that you can not know beforehand! Backgrounds are for removing really constant background signals.

Thresholds do not affect colocalization, but only the way colocalization maps are represented on the screen and objects splitted apart from each other based on the locally calculated colocalization. What affects colocalization, as explained in the colocalization theory at the SVI-wiki, are the backgrounds. In the computation of Manders coefficients the background values act like in this example: to the computation of  $m_1$  only pixels in R contribute when their corresponding pixel in G is above the background.

### Read more

More information is available at

<http://support.svi.nl/wiki/ColocalizationAnalyzer>

## 6. The PSF Distiller

Starting from version 2.5.7 Huygens Essential is optionally equipped with the *PSF Distiller* to allow you to measure your microscope's PSF. Measured PSFs allow you to improve deconvolution results and may also serve as a quality test for your microscope. It is a calibration of the microscope, in the sense of relating a physical known object with what the microscope actually measures.



The PSF Distiller is an optional component enabled by a `f` flag in the license string. Please refer to the *License string details* on page 49.

The PSF Distiller is an intelligent wizard based tool that is able to measure a PSF from one or more images of fluorescent beads, each containing one or more beads. It is also able to distill multi-channel PSFs from information collected from multi wavelength beads, or assemble a multi-channel PSF from single channel PSFs.

The PSF Distiller wizard leads you through four stages to carry out the following tasks:

- Preprocessing 'P'-stage to check microscopic image parameters. This stage is similar to the deconvolution 'P' stage.
- Stage 1: Averaging stage. In this stage all channels of the image are searched for beads that meet the selection criteria. After each successful or unsuccessful search there is the possibility to load in additional bead images, or go to the next stage.
- Stage 2: Distill stage. In this stage the PSF is measured from the averaged beads, for all available channels.
- Stage 3: Assembly stage. In case you want to combine results from earlier Distill runs with the current result to obtain a multi channel PSF you can add the earlier result here. It is also possible to add single or multi channel earlier results to a current multi channel result.




Figure 30. **The PSF Distiller.** After adding one or more bead images, an 'average bead' image is displayed in the Accumulator thumbnail. After pressing the DISTILL button, a PSF will be generated in a PSF thumbnail.

In the section *The Distiller stages* below, the stages will be discussed in detail.

### Beads for PSF distillation

The diameter of the beads should be in the order of the Half Intensity Width of the expected PSF. Larger beads will reduce the accuracy of the Distiller; smaller beads yield insufficient signal for accurate stacking in the averaging procedure, also resulting also in reduced accuracy. Typically beads with a diameter of 160 nm perform very well for many types of microscopy.



Beads should be recorded with **the same** microscopic parameters that you will use later to image your specimens. Please find more information about beads for PSF measurements on the SVI-wiki, at <http://support.svi.nl/wiki/RecordingBeads>.

## The Distiller stages

### Starting the Distiller

After launching Huygens Essential open the first bead image with **FILE > OPEN**. If your license includes the PSF Distiller option the **START PSF DISTILLER** button appears in the workspace; press it. One or more *Accumulator* images will be created into which later on the averaged beads will be kept. You now enter the 'P' stage.

### Parameter check stage

This stage is similar to the deconvolution 'P' stage and allows you to check all relevant microscopic parameters, in particular the sampling density.



Do not use undersampled bead images!

If any of the entry fields for the sampling density turn orange or worse, red, the data is unusable!

The parameters of bead images loaded at a later stage should match the ones you establish in this stage; if they do not you will be warned.

### Averaging stage

After setting the diameter and estimated Signal to Noise Ratio (SNR) of the beads, the image is searched for beads which meet all selection criteria:

- A bead should not be too close to another bead. If a bead is too close to another bead their images will interfere with each other. In widefield bead images this is quite problematic due to the large size of the blur cone. Fortunately, widefield PSF can be derived from single bead images.
- A bead should not be too close to an image edge. After all, another bead might be located just over the image edge.
- The intensity of a bead should not deviate too much from the median intensity of all beads. If it is brighter then it may be a cluster of two or more beads. If it is dimmer then it is not likely to be a bead. In both cases the object geometry is unknown so they are unusable.

If for some reason no beads are found, an explanation and some advice will be displayed in a popup. The software will try to find beads first with ideal selection criteria. If this does not yield a single bead then it will automatically retry with reduced inter bead distance criteria.

### Confocal and two photon bead images

Images from 160 nm beads should look like smooth fuzzy blobs with hardly visible noise. Use the default SNR settings. If available, it is a good idea to average 2–5 beads. Two photon bead images may look slightly noisy. If so, set the SNR to 20 and average 4–10 beads. To load additional bead images press **LOAD ADDITIONAL** after completing a search.

### Widefield bead images

Images from 160 nm beads should look like smooth fuzzy blobs with no visible noise. Use the default SNR settings. It is not necessary to average any more beads.

After you have added the last bead recording press `DISTILL` to enter the Distillation stage.

### Distillation stage

The distillation stage usually requires no user intervention, though in some cases a popup will be displayed with a question or a warning. All channels will be processed automatically.

### Assembly stage

In the assembly stage you can add a previously obtained PSF image file as a channel before or after the current result. Press `ADD CHANNEL` if you want to do this. The Distiller will compare the microscopic image parameters of the selected PSF image and check its content. In case there are problems the software will ask you to decide between ignoring the differences and discarding the selected file.

If you do not want to add channels press `NEXT` and then `DONE`. This leaves you with a loaded measured PSF for immediate use in a subsequent deconvolution run.



We advise you to measure the PSF of your system for each recording situation and certainly after each maintenance job in which the optics or scanning device was serviced.

Read more on the SVI-wiki, at <http://support.svi.nl/wiki/PsfDistiller>.

---

## 7. Establishing image parameters

The deconvolution algorithm needs to know some of the parameters describing your image acquisition. These are not too many, but you should carefully determine them. They are explained in this chapter.

---

### Image size

The amount of computing time involved in deconvolving images is more than proportional to the image size. It is therefore sensible to limit the data size as much as possible. With widefield images we recommend to not record planes below and above the object which only contain blur. Huygens Essential does not need these planes to restore your object. Since the blur in these planes might be affected by hard to correct bleaching they might even reduce the quality of the deconvolution result. In any case, never crop your objects of interest. As a rule of thumb, leave about one extra micron above and below your objects.



## Brick wise processing

Deconvolving images requires much computer memory because all computations are done in 32-bit floating point format, and because several extra (hidden) images are needed to store intermediate results. To reduce the memory requirements Huygens Essential will split your images into bricks, deconvolve the bricks sequentially, and fit the bricks together in a seamless fashion. Brick wise processing is an automatic feature of Huygens Essential. To find out the best number of bricks, let the software run in automatic mode for splitting. It will consider many options and go for the best one.

## Signal to Noise Ratio (SNR)

The Signal to Noise Ratio (SNR) can be estimated from the quality of the image. In Figure 31 you find some examples of recordings where different noise levels were added to an original (restored) image.




Figure 31. **Examples of different SNR values.** Same image, different SNR values. Top left: original image. Top right: image with SNR = 30. Bottom left: image with SNR = 15. Bottom right: image with SNR = 5.

If you want to calculate this SNR with more detail, you must remember that we **define the SNR** as the square root of the number of photons in the brightest part of the image.

Estimating the SNR in noisy images is surprisingly easy. It is based on the idea of establishing the voxel intensity value  $s$  corresponding to a single photon 'hit', by looking for such an event in dark areas of the image. Knowing the intensity value  $M$  of the brightest voxel in the image, one can now calculate how many photons are involved in it. Thus the SNR is

$$SNR = \sqrt{\frac{M}{s}}$$

If a significant blacklevel is present, it should be subtracted from these values before making the quotient. See the paragraph *Blacklevel* below for more details.

In low noise images it is much more difficult, as you no longer see single photon events so easily. Fortunately, in such cases the establishment of the precise SNR is not very important for the restoration method, and a rough estimation based on the appearance of the image is usually enough (see Figure 31).

Read more in <http://support.svi.nl/wiki/SignalToNoiseRatio>.

## Blacklevel

Figure 32 shows the histograms of three synthetic images. At left an image homogeneously filled with the value 5. At the middle we applied Poisson noise as if the




Figure 32. *Histogram of images with various blacklevel values.*

image was build by a CCD camera. At the image on the right a value of 20 was added to simulate electronic shift. This shift is called 'blacklevel'. A large blacklevel value will reduce the effective dynamic range of your microscope, but will do no harm to the deconvolution since it is automatically accounted for in the background estimation stage. However, it is also possible that the blacklevel is negative. In the image histogram this will show up as a spike on the left.

## Sampling densities

The **sampling density** is the number of recorded samples (voxels) per unit volume when acquiring the image. It is a microscopic parameter that describes the conditions of your image acquisition, established by the way you configure your microscope (usually by the *zoom* factor). The *ideal* sampling density depends on the system optics and is determined by the Nyquist rate: it is recommended that acquisition is done according to it as much as possible. If this was not achieved, the **actual** sampling distances must be used in any case as parameters when doing deconvolution. Imaging with one sampling density but using different values for deconvolution will produce wrong image restorations.

It is very important for the quality of a deconvolution result that all information generated by the optics of the microscope is captured in digital form. It can be shown that if the sampling density is higher than a certain value all information about the object is captured. We will call this value the critical sampling distance, corresponding to the **Nyquist rate**. Apart from practical problems like bleaching, acquisition time and data size there is no objection at all against using a smaller sampling distance than the critical distance, to the contrary.

Figure 33 shows the dependency of this critical sampling distance on the *NA* for specific wavelengths. To apply the plot of Figure 33 to different wavelengths you can simply scale the vertical axis with the wavelength. Example: you are working with a WF microscope with *NA* 1.3 at emission wavelength 570 nm. From the plot you read that the critical lateral Nyquist sampling distance at 500 nm emission is 95 nm, so in your case this becomes  $570 / 500 \times 95 \text{ nm} = 108 \text{ nm}$ .

In the confocal case it is the excitation wavelength which determines the Nyquist sample distance. In theory the pinhole plays no role, but larger pinholes strongly attenuate fine structures at the resolution limit. Therefore, as a rule of thumb, with a common pinhole diameter of 1 Airy disk the lateral critical sampling distance may be increased by 50% with negligible loss of information. In cases where the pinhole is much larger, the lateral imaging properties much resemble those of a WF system and the sampling distance can be set accordingly. We do not recommend to increase the axial sampling distance appreciably beyond the critical distance

In a multi-photon excitation microscope, it is the excitation wavelength divided by the photon count which determines the sampling.

Read more on the SVI-wiki: <http://support.svi.nl/wiki/NyquistRate>.




Figure 33. **Critical sampling distance vs. NA.** The curves above show the critical sampling distance in axial and lateral directions for wide-field and confocal microscopes. The emission wavelength in both cases is 500 nm.

## Computing the backprojected pinhole radius



Throughout the Huygens Professional and Huygens Essential pinhole sizes of confocal systems are specified as the **backprojected radius in nanometers**  $r_b$ . 'Backprojected' means the size of the pinhole as it appears in the specimen plane: the physical pinhole radius  $r_{phys}$  divided by the total magnification of the detection system. This total magnification is the product of the (variable) objective magnification times a fixed internal magnification:

$$r_b = \frac{r_{phys}}{m_{obj} m_{system}} \quad (1)$$

where  $m_{obj}$  is the magnification factor of the objective and  $m_{system}$  is the fixed magnification of the system.

The equations that you can find in the next pages for different microscopes are intended to orientate the user in finding out the backprojected value, but the idea is always the same: given the *diameter* of the real pinhole  $d$ , we apply a factor for unit conversion (to obtain the *radius* in nm), and some dividing numbers that take account of the magnification of the microscope. These include both the objective and the intrinsic system magnification. In some microscopes with pinhole shapes other than circular, a geometrical correction will also be needed.

### Airy disk as unit for the backprojected pinhole

Some confocal microscopes report the pinhole size (diameter) with the Airy disk (diameter) as unit. The backprojected pinhole radius can then be computed with:

$$r_b = \frac{0.61 \lambda_{ex} N_{Airydisks}}{NA} \quad (2)$$

with  $NA$  the numerical aperture of the lens,  $N_{Airydisks}$  the number of Airy disks, and  $\lambda_{ex}$  the excitation wavelength. In principle using  $\lambda_{ex}$  is not correct because the Airy diffraction pattern is formed by the emitted light. However, we suspect microscope manufacturers prefer to use the excitation wavelength because it is better defined and does not depend on settings of devices like adjustable band filters. For this reason in the formula above we too use  $\lambda_{ex}$ .

Note that this relation bypasses the need to know internal system and lens magnifications.

### Converting from integer parameter


Unfortunately, quite a few commercial microscopes do not report the physical pinhole size or the Airy disk size. Instead, often an integer size parameter is specified with a range 0–255 (8-bit). Matters are further complicated by the use of non-circular pinholes. To compensate this we introduce a shape factor  $c_{shape}$  which takes care of the conversion from *size* (diameter or edge size) to *radius*. The following formula can be used to translate the 8-bit machine number into the backprojected pinhole radius:

$$r_b = 1000 \ c_{shape} \frac{(P_8/255) (s_{max} - s_{min}) + s_{min}}{m_{obj} \ m_{system}} \quad (3)$$

where  $P_8$  is the 8-bit machine number for the pinhole,  $s_{min}$  is the maximal pinhole size in *micrometers*,  $s_{min}$  is the minimal pinhole size in micrometers. The factor 1000 is to convert from microns to nanometers.

### Shape factor

The shape correction from a square to a circular pinhole is based on setting the area of the square pinhole equal to the area of the replacing circular pinhole.



$$d^2 = \pi r^2 \text{ and thus } r = \frac{d}{\sqrt{\pi}} \quad (4)$$

For a square pinhole  $c_{shape} = 1/\sqrt{\pi} = 0.564$ . For circular pinholes  $c_{shape} = 0.5$  to just convert from diameter to radius.

### Airy disk as unit for the backprojected radius of a square pinhole

The relation between the edge size of square pinhole in Airy disk units and the backprojected radius is a combination of equations 2 and 4. Taking into account that equation 2 already converts a diameter into a radius, we get:

$$r_b = \frac{0.69 \lambda_{ex} N_{Airydisks}}{NA} \quad (5)$$

with  $N_{Airydisks}$  the number of Airy disks.

### Computing the backprojected pinhole distance in Nipkow spinning disks

As is the case for the backprojected pinhole diameter, the distances between the pinholes in spinning disks must be divided by the system magnification. For the most used Yokogawa, for example, the pinhole physical distance is 250  $\mu\text{m}$ . You can check this by stopping the disk. So with an  $\times 100$  lens the backprojected distance is about 2.5  $\mu\text{m}$ . If an extra zoom lens is placed between the disk and the sample, its magnification must be also considered.

---

## Pinhole radius tables

In the case of the Leica TCS, the Biorad MRC, the Olympus Fluoview and the Zeiss LSM microscopes the pinhole geometry and system magnification is known resulting in the conversion formulas from the following tables:

### Leica confocal microscopes

#### TCS 4d, SP1, NT

The size of the square pinhole is given as an 8-bit number which maps to the physical pinhole radius given in the following table:

TCS4d	Range begin	Range end
Reported parameter ( $P_8$ )	0	255
Diameter (micron)	20	630 (earlier reported as 500)
Pinhole geometry	square	
System magnification	4.5	
Backprojected pinhole radius (nm)	$\frac{295 P_8 + 2508}{m_{obj}} \quad (\text{Eq. 6})$	

Table 1. Leica TCS4d pinhole parameters. This table is believed to be valid for the SP1 and NT models too.

If the pinhole is specified in Airy disk units, see *Airy disk as unit for the backprojected radius of a square pinhole* on page 43 for information on how to convert to a backprojected radius.

### TCS-SP2

The Leica TCS-SP2 has a system magnification of 3.6. However, the size of its square pinhole is usually specified in Airy disk units making it independent of the actual overall magnification. It is dependent though on wavelength and  $NA$  of the objective. See *Airy disk as unit for the backprojected radius of a square pinhole* on page 43 for information on how to convert to a backprojected radius.

## Zeiss confocal microscopes

Zeiss LSM410 inverted	Range begin	Range end
Reported parameter ( $P_8$ )	0	255
Diameter (micron)	0	1000
Pinhole geometry	square	
System magnification ( $m_{sys}$ )	2.23	
Backprojected pinhole radius (nm)	$\frac{992 P_8}{m_{obj}} \quad (\text{Eq. 7})$	

Table 2. Zeiss LSM410 (inverted) pinhole parameters.

Zeiss LSM510	
Diameter (micron)	$d$
Pinhole geometry	circular
System magnification ( $m_{sys}$ )	3.33
Backprojected pinhole radius (nm)	$\frac{10^3 d}{2 m_{sys} m_{obj}} \quad (\text{Eq. 8})$

Table 3. Zeiss LSM510 pinhole parameters.

## Olympus confocal microscopes

Olympus Fluoview					
Reported parameter	1	2	3	4	5
Diameter $d$ (micron)	0	100	150	200	300
Pinhole geometry	square				
System magnification ( $m_{\text{sys}}$ )	3.426 / 3.8 (FV500)				
Backprojected pinhole radius (nm)	$\frac{10^3 d}{2 \sqrt{\pi} m_{\text{sys}} m_{\text{obj}}} \quad (\text{Eq. 9})$				

Table 4. Olympus Fluoview pinhole parameters.

## Biorad confocal microscopes

Biorad	
Reported parameter ( $P_g$ )	n.a.
Diameter $d$ (mm)	0–8
Pinhole geometry	circular
System magnification ( $m_{\text{sys}}$ )	53–83, reported 60 for the Radiance, 53 for the 1024.
Backprojected pinhole radius (nm)	$\frac{10^6 d}{2 m_{\text{system}} m_{\text{obj}}} \quad (\text{Eq. 10})$

Table 5. Biorad MRC600/1024/Radiance pinhole parameters.

### Checking the Biorad system magnification

The Biorad MRC500/600/1024 microscopes have a very high magnification in the detection system. The fixed system magnification is, according to Pawley<sup>9</sup>,  $53 \times m_{\text{tube}}$ , with  $m_{\text{tube}}$  between 1.0 and 1.56 (a factor 1.25 for the fluorescence attachment and also a factor 1.25 for the 'DIC' attachment). The factor of 53 includes the 8× "eyepiece" just below the scan head, but does not include that variability in magnification due to the variations in tube-length that are result from the aligning the system.

The high system magnification allows you to view the diffraction pattern (Airy disk) at the pinhole plane directly by eye. To enable you to verify the correctness of equation 10 in the table above for your instrument we outline the way the system magnification was derived:

In a Biorad MRC600 with a 1.3 60× objective, the Airy disk has a diameter of around 2–2.5 mm at the pinhole plane. The diameter of the first Airy zero ring is 7.6 lateral optical units (o.u.), using the following transformation to express a distance  $x$  in dimensionless o.u.,

$$v = x \frac{2 \pi}{\lambda} n \sin \alpha \quad (11)$$

where  $n \sin \alpha$  is the numerical aperture  $NA$ . In the system above an o.u. is  $0.3 \pm 0.033$  mm. At the specimen plane (backprojected) a lateral o.u. is in this case around 61 nm. The total magnification

<sup>9</sup> Pawley, J. B. Handbook of biological confocal microscopy, 2<sup>nd</sup> ed. 1995. Plenum Press, New York and London. ISBN 0-306-448262. Page 30.



appears thus to be 4918, the system magnification  $82 \pm 9$ . This value corresponds well with the largest possible system magnification for the MRC600.

### A supplied calibration curve

If a calibration curve was supplied with your microscope best use that curve to convert the displayed setting to a physical size and from there convert to the backprojected radius.

---

## 8. Improving the quality of your images

In this chapter you can find some basic suggestions on how to acquire better microscope images. They are based on frequent problems that we find on users data. These recommendations go on the direction of obtaining the greatest quality images from your microscope, from the point of view of acquiring as much information as clean as possible. This alone worths the effort, but it will also be very valuable for the deconvolution afterwards.

Some basic guidelines to improve the deconvolution results are also included. You can find more information at the SVI-wiki, at <http://support.svi.nl/wiki>.

---

### Data acquisition pitfalls

Some of the most common problems arising when acquiring images are addressed in this section.

#### Refractive index mismatch

A mismatch between the refractive index of the lens immersion medium and specimen embedding medium will cause several serious problems:

- Geometrical distortion: the fishtank effect  
Objects will appear elongated in the microscope. Huygens Essential takes this into account when calculating the theoretical PSF, but it will not modify the image geometry.
- Spherical aberration (S.A.)  
S.A. will cause the oblique rays to be focused in a different location than the central rays. The distance in this focal shift is dependent on the depth of the focus in the specimen. If the mismatch is large, e.g. when going from oil immersion into a watery medium, the PSF will become asymmetric at depths of already a few micron. Especially harmful for WF deconvolution. Workaround: keep the *z*-range of the data as small as possible. Solution: use a water immersion lens.
- Total internal reflection  
When the lens numerical aperture *NA* is larger than the medium refractive index *n* total internal reflection will occur, causing excitation light to be bounced back into the lens and limiting the effective *NA*.

If S.A. is unavoidable, you can still improve the image during restoration using an adaptive point spread function. See *Spherical aberration correction* on page 48.

#### Clipping

The light intensities from the microscopic object are converted to electrical signals that pass an adjustable amplifier. Also an electrical DC component can be added or subtracted by the microscope operator. The electrical signal may thus range from negative to highly positive. These electrical signals must be converted to numbers processed by the computer. This converting stage is

done in the CCD camera and its electronics. Most CCD cameras have an 12-bit converter limiting the output numbers to a range of 0 to 4095. Negative input signals are usually converted to 0 while positive input values exceeding some value are all converted to 4095 (clipping): information in the clipped samples is *lost*.

In practice: be suspicious if you find in your image voxel intensity values at the extremes of the numerical range, probably clipping occurred.

Read more in <http://support.svi.nl/wiki/ClippedImages>.

## Undersampling

One of the rules of measurement that is often overlooked is that one takes too few *xy* slices from the microscopic object. In that case the sampling distance is too large (too few samples: under-sampled) which leaves you with a 3-D stack with hardly any relation between the adjacent planes. It is important to know how the sampling conditions should be established in order to recover an image from the sampled values. How you should sample your object depends on your microscope type -WF or confocal- and on the microscope parameters used, like the numerical aperture and wavelength of the light. The correct sampling distances can be estimated as explained on *Sampling densities* (page 40) or calculated using the formulas as given in *The Nyquist rate* from the **Huygens Deconvolution Recipes**<sup>10</sup>.

## Do not undersample to limit photodamage

Some times undersampling is done to limit photodamage to live cells. If photodamage plays a role it is better to distribute the available photons over more pixels, resulting in an apparently noisier image, than putting the photons in fewer pixels to get a low noise, but undersampled, image. Of course there are limits, but a fair trade-off can be often found.

It is better to record 10 separate noisy slices 100 nm apart than slices on 1000 nm each averaged 5 times in order to reduce noise. See also *A typical example* in *Sampling an image* from the **Huygens Deconvolution Recipes** or in the above mentioned Internet link.

## Bleaching

Bleaching is a practically unavoidable phenomenon in fluorescence microscopy. Because the image planes are acquired sequentially, bleaching will vary along the *z* direction. Assuming it is not strong it will not affect deconvolution results on confocal or two photon images. But in WF deconvolution bleaching is more of a problem. Fortunately, usually the bleaching in WF images can be corrected quite easily. Huygens will do so automatically. However, if the bleaching is strong the correction might not be perfect, resulting in lower quality deconvolution results.

## Illumination instability

Some WF systems are equipped with unstable arc lamps. Huygens will correct this instability, but when the instability is severe it cannot do so sufficiently.

## Mechanical instability

Mechanical instability can take many shapes, for example:

- *Vibrations sometimes seen in confocal images*. They may seriously hamper deconvolution
- *z-stage moves irregular or with sudden jumps*. Fatal for confocal or WF deconvolution.
- *Specimen moves*. If in WF data the object can clearly be seen moving when slicing along over a few micron in *z* this will cause problems for the deconvolution. Best cause of ac-

<sup>10</sup> Also available on-line in <http://support.svi.nl/wiki/NyquistRate>

tion, apart from speeding up acquisition, is limiting the z-range of the data as much as possible. Confocal data of moving specimen causes less problems.

### Thermal effects

Thermal effects are known to affect calibration of the z-stage, especially if piezo actuators without feedback control are used. In particular harmful for WF data.

### Internal reflection

At high  $NA$  the angle of incidence of the most oblique rays can be close to 70 degrees. When a ray has to cross the cover-glass to medium interface at such an angle total reflection may occur. To be precise, total reflection occurs when the  $NA$  of your lens is higher than the refractive index of the embedding medium. This will reduce the effective  $NA$  of the lens.

---

## Deconvolution improvements

### Acquire an experimental PSF

A point spread function (PSF) is the image of a single point object. The degree of spreading (blurring) in the image of this point object is a measure for the quality of an optical system, and the imaging in a fluorescent microscope is completely described by its PSF. Although in many cases a theoretically calculated PSF very well matches the real one, ideal theoretical calculations can not predict actual misalignments or other problems inside the optical path. Therefore it is always recommendable to measure an experimental PSF and, if it is very different from the ideal one, use it for deconvolution instead of the theoretical one.

Obtaining an experimental PSF is based on recording the image of a well know probe, as close as the ideal subresolution light source as possible, and distilling a PSF from it using *The PSF Distiller* as explained on page 36.

### Spherical aberration correction

The Huygens software automatically adapts the **theoretical** point spread function (PSF) to the sample depth to correct for spherical aberration, in case of a refractive index mismatch. For that, the program considers that the coverslip is placed at the bottom of the image, at lower Z coordinates (as in an inverted microscope). If your image suffers from a refractive index mismatch you should adapt the image to this condition before the restoration: the Huygens Essential has a MIRROR ALONG Z tool, to flip the image in case your coverslip is on the top.

In case your image suffers from high spherical aberration, it might be better to use a theoretical PSF with this depth-dependent correction than an experimental one.

Read more in <http://support.svi.nl/wiki/SphericalAberration>.

### Improve the estimated parameters

This is an usual procedure to improve the estimated parameters for deconvolution: you start out with a signal-to-noise ratio (SNR) derived either from an estimate of the number of photons in the image or simply from a visual inspection (as explained on page 39). Subsequently you do a restoration run and inspect the result for artifacts and residue background. If you are confident all is fine you resume the restoration with a higher SNR setting (say 30%-50%) and perhaps a higher background. The software will ask you to continue were you left off (keeping improving the image, quite recommendable) or to start from the raw image again with the new parameters. A new result

will be generated to compare with the previous one, for instance using the Twin Slicer. You can repeat this several times, and at the end you will be asked to select the best result as the final one.

If you have done this a previous time for a similar image then of course you just use the values established then!

Find a detailed example in the following wiki article:

<http://support.svi.nl/wiki/DeconvolutionProcedure>.

## 9. Appendix

### License string details

You can have detailed information about your license strings by going to **HELP > LICENSE**, select the license string you want to inspect, and click **EXPLAIN LICENSE**.

This is a License String example:

```
HuEss-2.7-wcnp4-d-ft-demnps-2008Sep30-33dfa680a8402167-
{info@svi.nl}-c338a4da57cb342151d3
```

As from Huygens Compute engine version 2.4 the license strings have the following format:

```
<prod>-<vers>-<micr>-<server>-<flags>-<lock>-<date>-<sysid>-
<email>-<checksum>
```

Where:

<prod> can be one of HuEss, HuScript, HuPro, FluVR.

<vers> is the version number which is 2.4 or higher.

<micr> consists of one or more characters representing microscope types:

w: Widefield, enables you to deconvolve WF data

c: Confocal, enables you to deconvolve confocal data

n: Nipkow Disk, enables you to deconvolve data from Nipkow (spinning) Disks

p: multi-photon, for 2-(or more) photon microscopes

4: 4-pi, for experimental microscopes

<server> can be one of:

d: desktop (1 or 2 CPU's)

s: small server (3–8 CPU's)

m: medium server (9–32 CPU's)

l: large server (33–64 CPU's)

x: extreme server (65 or more CPU's)



Remark for Linux and Mac OS X multi-processing users:

d: desktop (1 thread)

m: medium server (max 2 threads)

l: large server (max 4 threads)  
x: extreme server (5 threads and more)

<extra flags> which enables additional functionality:

f: Distiller option, enables obtaining PSF's from measured beads.  
t: Time series option, enables deconvolution of time series, either 2D-time or 3D-time  
v: Surface Renderer visualization option.  
C: Colocalization Analyzer option.  
A: Object Analyzer option.

<lock> Set of license policy or locking flags having of one or more characters:

d: License expiry date  
e: Email address  
m: Hardware system ID  
n: Number of processors  
p: Processor type  
s: Processor details

<date> is the end of license or end of support date in the format eomYYYYMMDD

<sysid> is 16 hexadecimal digits of system identification

<email> is the email address of the customer who bought this license in the format {email address}

<checksum> is 20 hexadecimal digits of license checksum

When the d flag (license expiry date) is specified, then after the locking flags comes the expiry date. Without this flag the expiry date will be preceded with the letters eom, which turns the date into an end date of the current support contract. Currently the e flag is always present. If it is after the system ID an email address surrounded by braces is given. This is the email address of either the creator of the license, or the sales person in case of a temporary license, or the customer to whom the license was sold. This is used for informational purposes only. It does not limit the license in any way. If the m flag is present then the hardware system identification ID must match the one from the running system. Similar for the n flag. It requires that the number of processors as specified by the system ID matches the number in the current system. The p requires a match on processor type. For PCs this consists of the name of the processor vendor like Intel or AMD, the processor family like Pentium, and the processor model. For Irix machines this is the IP level of the processor board which can be found in the first line of the output of the hinv command. The s flag is only meaningful for PCs. It enables additional matching on the processor core 'stepping number' which is a notation from Intel (like a revision number) that indicates the changes or improvements inside the CPU's instructions set.

---

## Questions

### Where can I find support on the web?

See *Addresses and URLs* on page 53.

### What does the quality factor mean while running Huygens?

Deconvolution as it is done in Huygens Essential hinges around the idea of finding an as good as possible estimate of the object that is imaged by the microscope. To assess the quality of an estimate, Huygens Essential computes the image of each estimate as it would appear in the microscope

and compares it with the measured image. From the difference a quality factor is computed. The difference is also used to compute a correction factor to modify the estimate in such a way that the corrected estimate will yield a better quality factor. The quality factor as reported by Huygens Essential is a measure relative to the first estimate and therefore a number greater or equal to 1. If the increase in quality drops below a threshold the iterations are stopped. See also *Quality threshold* on page 15.

### Can I deconvolve a TIFF series?

Yes, if the series is a numbered series like: `slice001.tif`, `slice002.tif`, ..., `slice0nn.tif` Huygens Essential will read the series into a single 3D image. Because TIFF's usually carry no additional microscopic information, check the parameters carefully.

### TIFF file series naming convention

If you have TIFF images to be read into the Huygens Professional or the Huygens Essential you should know about the naming convention used.

If you select a file from a numbered series, the selected file and the following files will be interpreted as x-y planes of a 3D stack and read into a 3D image of suitable size and channel configuration.

A one-channel 3D images only go with numbers:

As an example a dataset called `c` with 32 slices numbers as follows:

```
c000.tif
c001.tif
...
...
c031.tif
```

If you wish to work on the complete `c` stack you only have to select `c000.tif` while opening. If you select file `c020.tif` the first 20 slices will be skipped.

Numbered series without the TIFF extension like

```
c04
c05
...
...
c18
```

are not read in as a series.

Huygens Essential and Professional read and write TIFF series with Leica style numbering if there are more channels (different wavelengths), slices or frames(in time).

An image of four slices and two frames is named with Leica style numbering as follows:

```
c_t00_z000.tif
c_t00_z001.tif
c_t00_z002.tif
c_t00_z003.tif
c_t01_z000.tif
c_t01_z001.tif
c_t01_z002.tif
c_t01_z003.tif
```

And an image sTCh of four slices, three frames and two channels:

```
sTCh_t00_z000_ch00.tif
sTCh_t00_z000_ch01.tif
sTCh_t00_z001_ch00.tif
sTCh_t00_z001_ch01.tif
sTCh_t00_z002_ch00.tif
sTCh_t00_z002_ch01.tif
sTCh_t00_z003_ch00.tif
sTCh_t00_z003_ch01.tif
sTCh_t01_z000_ch00.tif
sTCh_t01_z000_ch01.tif
sTCh_t01_z001_ch00.tif
sTCh_t01_z001_ch01.tif
sTCh_t01_z002_ch00.tif
sTCh_t01_z002_ch01.tif
sTCh_t01_z003_ch00.tif
sTCh_t01_z003_ch01.tif
sTCh_t02_z000_ch00.tif
sTCh_t02_z000_ch01.tif
sTCh_t02_z001_ch00.tif
sTCh_t02_z001_ch01.tif
sTCh_t02_z002_ch00.tif
sTCh_t02_z002_ch01.tif
sTCh_t02_z003_ch00.tif
sTCh_t02_z003_ch01.tif
```

### **Can I deconvolve a single plane widefield image?**

Yes. Single plane WF deconvolution works because the data is extrapolated into a region above and below the plane spanning typically between 10–20 planes of 100–300 nm sampling in *z*. The software generates an appropriate PSF.

### **Can I deconvolve a single TIFF image?**

Yes. Huygens Essential treats the image as the only known plane of a 3D stack and proceeds as usual. Set the *z*-sampling distance to the Nyquist value: see *Sampling densities* on page 40.

### **How do I generate a debug log?**

See <http://support.svi.nl/wiki/DebugMode>.



---

## Addresses and URLs

### Where can we be reached?

#### Scientific Volume Imaging b.v.

Alexanderlaan 14  
1213XS Hilversum  
The Netherlands

You can call us directly by phone:

+31 35 6859405 or ++653 345445,  
or fax us at: +31 35 6837971,  
or email us at: [info@svi.nl](mailto:info@svi.nl)

<http://www.svi.nl>

#### Distributors

A list of distributors can be found on our web site:

<http://www.svi.nl/company/distributors.php>

### Support and FAQ

#### Knowledge database (FAQ)

It is good to know that an extensive support base is available where you may find answers to questions that come to front while reading this document. FAQ's are available for many items such as deconvolution and general microscopy, installation, memory management, visualization, file formats, platforms and reported bugs:

<http://support.svi.nl>

On this web site you will also find a form to submit questions to SVI's support team.

#### The SVI-wiki

The SVI-wiki is a rapidly expanding public knowledge resource on 3D microscopy and deconvolution. Based on the WikiWikiWeb principle, it is open to contributions from every visitor. In addition it serves as a support medium for SVI customers and relations to discuss different aspects of SVI's Huygens Software. The SVI-wiki can be found in

<http://support.svi.nl/wiki>

#### Starting points

Here you can find a selection of good starting points in the SVI-wiki to learn more about the Huygens Software and microscopical imaging:

- On the parameters describing your **imaging conditions** (sampling, numerical aperture, pinholes...) see <http://support.svi.nl/wiki/MicroscopicParameters>

- On the **restoration parameters** (signal to noise ratio, backgrounds, quality criteria...) used by the deconvolution algorithms, see <http://support.svi.nl/wiki/RestorationParameters>
- A step by step example on how to **tune these parameters** to achieve the desired restoration results: <http://support.svi.nl/wiki/DeconvolutionProcedure>
- Very **important factors** (sampling, clipping...) on image acquisition and restoration: <http://support.svi.nl/wiki/ImportantFactors>
- Typical **acquisition pitfalls** (spherical aberration, undersampling, bleaching...) explained: <http://support.svi.nl/wiki/AcquisitionPitfalls>
- On recording beads to **measure a PSF**: <http://support.svi.nl/wiki/RecordingBeads>
- **Tutorials** and detailed information on using the different software **tools** (restoration, visualization, analysis, programming...) can be found in <http://support.svi.nl/wiki/Tutorials>
- If you are **sending your images** to SVI to get further assistance, please follow the instructions detailed in <http://support.svi.nl/wiki/SendImagesToSvi>.
- Do you want to **edit the wiki**, contributing to enhance its contents? Read <http://support.svi.nl/wiki/HowToUseWiki>
- With the **search tools** you can find much more information on the topic you want: <http://support.svi.nl/wiki/FindPage>

## Quick reference

Write here your support data in order to have it always available:

<b>SVI support number</b>	
<b>Login for the wiki and the image database</b> (do not write your secret password!)	

For further information see <http://support.svi.nl/wiki/SupportNumber> and <http://support.svi.nl/wiki/WhyLogin>.

# Alphabetical Index

## A

acquisition pitfalls 46  
acquisition time 40  
Airy disk 41-45  
Analyzer  
    Colocalization Analyzer 32  
    Object Analyzer 28  
anchor 32  
animation 23, 26  
artifacts 46  
average background 13

## B

background 34, 35  
backprojected distance 9, 43  
backprojected pinhole 9, 41-45  
batch processing 18  
bead 7, 14, 36-38, 50  
Biorad 16, 43, 45  
blacklevel 12, 40  
bleaching 15, 38, 40, 47  
brick 14, 39  
brightfield images 8

## C

caliper 31  
channel removal 12  
circular pinhole 42, 43  
clipping 12, 46, 47  
CMLE 19  
Colocalization Analyzer 28, 32, 50  
colocalization coefficients 34  
colocalization map 34  
color 27  
color mode 21  
color of the channels 9  
compress color mode 22  
compute engine 1, 49  
computing time 38  
confocal 37, 41, 49  
contrast 8, 21, 22  
correction 15, 42, 47, 51  
correction for z-drift 16  
critical sampling 9  
critical sampling distance 40, 41  
cropping 7, 11, 12, 31  
cyan background 9

## D

debug 52  
deconvolving 8, 38, 39  
demo image 7, 12  
distance 9  
distributors 53  
download 1, 4  
drift corrector 16

## E

email SVI 53  
emission 24, 25, 40, 41  
excitation 24, 25, 41, 42, 46

## F

FAQ 53  
Fast Classic Maximum Likelihood Estimation 14  
faxing SVI 53  
fishtank effect 46  
Fluorescence Volume Renderer 24  
FluVR 24, 49  
Frequently Asked Questions (FAQ) 53

## G

geometrical distortion 46

## H

histogram 7, 12, 13, 40  
hue 27  
HuEss 2, 49  
HuPro 49  
HuScript 49

## I

ICS file 16  
ICS2 file 16  
illumination instability 47  
image statistics 9  
Imaris-Classic file 16  
internal reflection 46, 48  
inversion of contrast 8  
isosurface 28

## L

labeling 28  
Leica 43, 44, 51  
Leica style numbering 51  
license string 2-4, 26, 28, 33, 36, 49  
linear color mode 22

## M

Maximum Intensity Projections 11  
Maximum Likelihood Estimation 14, 15  
mechanical instability 47  
MIP 11, 23, 26  
MIP Renderer 23  
MLE 14-16  
movie rendering 23, 26  
multi-channel 7, 12, 16, 17, 24, 36  
multi-photon excitation 37, 41, 47, 49

## N

NA 9, 40, 42, 44-46, 48  
Nipkow 9, 43, 49  
non-circular 42  
numbering of TIFF series 51  
numerical aperture 9, 42, 45, 47  
Nyquist rate 40, 47

## O

object 28  
object analysis 33  
Object Analyzer 28  
Olympus 45  
optical units 45  
orange background 9

## P

photon 37, 39, 47  
photon count 9  
pinhole 9, 41-45  
Point Spread Function (PSF) 14  
PSF 8, 14, 15, 36, 46, 52  
PSF Distiller 4, 7, 14, 36

## Q

QMLE 19  
quality threshold 15

## R

red background 9  
refractive index mismatch 46  
region of interest 29, 30  
Renderer  
    MIP 23  
    SFP 24  
    Surface 26  
rendering 23, 25-27  
restoration algorithm 19

ROI 29, 30

## S

sampling 13, 37, 40, 47  
saturation 12  
saving 16, 17  
Scientific Volume Imaging b.v. 53  
seed and threshold 30  
segmentation 28  
SFP renderer 24  
shape correction 42  
shape factor 42  
Signal to noise ratio (SNR) 15, 39  
Simulated Fluorescence Process 24  
SNR 15, 39  
specimen moves 47  
spherical aberration 14, 46, 48  
square pinhole 42-44  
statistics 9  
StuffIt Expander 1  
Surface Renderer 26, 50

## T

task report 17  
template 10  
thermal effects 48  
threshold 30, 35  
TIFF file 16, 27, 51  
TIFF series 26, 51  
time series 8, 15, 16, 22, 24, 26, 50  
total internal reflection 46  
Transparency 25, 27  
Twin Slicer 14, 16, 21, 22  
two photon 37, 47, 49

## U

undersampling 47  
unstable arc lamps 47  
updates 3

## V

vibrations 47  
Viewer  
    MIP Renderer 23  
    SFP Renderer 24  
    Surface Renderer 26  
    Twin Slicer 21  
violet background 9

## W

wavelength 9, 24, 25, 36, 40-42, 44, 47, 51  
website 4, 24  
widefield 38, 49

## Z

z-drift corrector 16  
z-stage moving 47  
Zeiss 43, 44

## .

.exe file 1  
.gif file 24  
.ics file 17  
.ids file 17  
.pic file 16  
.pkg file 1  
.rpm file 2  
.tar file 2  
.tardist file 2  
.xml file 11

## 2

2D images 8

## 4

4-pi 49